

Public Health Inspector's Guide to Environmental Microbiology Laboratory Testing



Current as of August 8, 2021

Public Health Ontario

Public Health Ontario is a Crown corporation dedicated to protecting and promoting the health of all Ontarians and reducing inequities in health. Public Health Ontario links public health practitioners, frontline health workers and researchers to the best scientific intelligence and knowledge from around the world.

Public Health Ontario provides expert scientific and technical support to government, local public health units and health care providers relating to the following:

- communicable and infectious diseases
- infection prevention and control
- environmental and occupational health
- emergency preparedness
- health promotion, chronic disease and injury prevention
- public health laboratory services

Public Health Ontario's work also includes surveillance, epidemiology, research, professional development and knowledge services. For more information, visit <u>publichealthontario.ca</u>.

How to cite this document:

Ontario Agency for Health Protection and Promotion (Public Health Ontario). Public health inspector's guide to environmental microbiology laboratory testing. Evergreen ed. Toronto, ON: Queen's Printer for Ontario; 2021.

©Queen's Printer for Ontario, 2021

Public Health Ontario is an agency of the Government of Ontario.



Authors

Vanessa G. Allen, MD, MPH Chief, Medical Microbiology Public Health Ontario

Sandra Edelsward, B.Sc., MLT Program Coordinator, Drinking Water Testing Public Health Ontario

Allana Murphy, B.Sc.
Senior Laboratory Lead - Environmental
Microbiology
Public Health Ontario

Anna Majury, DVM, PhD Clinical Microbiologist Public Health Ontario

Anne Maki, MLT Senior Manager, Clinical and Environmental Bacteriology Public Health Ontario

Acknowledgements

The authors wish to express their sincere appreciation for the effort and dedication demonstrated by Ontario's Public Health Units (PHUs) throughout the enhancements of the Public Health Inspector's Guide to Environmental Microbiology Laboratory Testing (the Guide). We would like to thank our provincial and federal partners for their feedback:

- Association of Supervisors of Public Health Inspectors of Ontario (Canada) (ASPHIO)
- Canadian Food Inspection Agency (CFIA)
- Canadian Institute of Public Health Inspectors (CIPHI)
- Ministry of the Environment, Conservation and Parks (MECP)
- Ontario Ministry of Agriculture, Food and Rural Affairs (OMAFRA)
- Ministry of Health, Environmental Health Policy and Programs, Health Protection and Surveillance Policy and Programs Branch (MOH-EH)

We also thank our colleagues at Public Health Ontario for their collaboration and support in the development of the Guide.

Disclaimer

This document was developed by Public Health Ontario (PHO). PHO provides scientific and technical advice to Ontario's government, public health organizations and health care providers. PHO's work is guided by the current best available evidence at the time of publication.

This document is intended to assist public health inspectors in decision-making by describing PHO Laboratory (PHOL) testing services for environmental microbiology. This document should not be considered inclusive of all proper methods of environmental testing or exclusive of other methods of environmental testing reasonably directed at obtaining the same results. The application and use of this document is the responsibility of the user. PHO assumes no liability resulting from any such application or use.

This document may be reproduced without permission for non-commercial purposes only and provided that appropriate credit is given to PHO. No changes and/or modifications may be made to this document without express permission written from PHO.

Note:

Regular updates will be made to the Guide and changes will be tracked in the revision history. It is the responsibility of the users to determine whether there have been any such updates or amendments, and PHO shall have no responsibility for providing notice of any such updates or amendments.

The electronic version of this Guide includes links to third party resources ("Other Resources") which are not maintained by PHO. PHO does not review, endorse or approve the Other Resources and is not responsible for its content. PHO will not be liable for any loss or damages arising from your use of the Other Resources."

Contents

ln	troduction	1
	Environmental Matrices Analysed by the PHO Laboratory	2
	Investigations and Analyses Not Available at the PHO Laboratory	3
	Other Resources	3
	Legal Considerations	4
	Legal Samples – Collection Instructions	4
	Revision History	5
Fo	ood Sample Analysis	6
	Food - Microbiological Analyses Test Menu	6
	Food Surveys or Studies	8
	Sample Selection	8
	Containers for Sampling	9
	Sample Preparation	9
	Sample Collection	10
	Sample Collection – Large Lots of Food	11
	Sample Collection – Small Quantities of Food (<200 grams)	11
	Completion of the Food Bacteriology Requisition	11
	Shipping Instructions	12
	Food – Laboratory Acceptance Criteria	12
	Food – Reporting and Interpretation of Environmental Microbiological Test Results	13
	Revision History	17
Eı	nvironmental Swabs	18
	Swab – Microbiological Analysis Test Menu	18
	Sample Selection	18
	Containers for Sampling	18
	Sample Preparation	19
	Sample Collection	19
	Environmental Swab – Laboratory Acceptance Criteria	20
	Environmental Swab – Reporting and Interpretation of Environmental Microbiological Test Results	21
	Revision History	.22

Water Sample Analysis	23
Water – Microbiological Analysis Test Menu	23
Important Information Related to Drinking Water Analysis	24
Chain of Custody	26
Water – Laboratory Acceptance Criteria	28
Water – Reporting and Acceptable Limits of Environmental Microbiological Test Results	30
Water Status Codes and Associated Interpretation Statements for Drinking Water Reports	34
Revision History	35
Legionella Investigations	36
Legionella - Microbiological Analysis Test Menu	36
Sample Selection	37
Containers for Sampling	39
Sample Preparation	39
Sample Collection - Swabs	40
Sample Collection - Water	40
Legionella – Laboratory Acceptance Criteria	41
Legionella – Reporting and Interpretation of Environmental Microbiological Test Results	42
Revision History	42
Glossary	43
Acronyms	51
Revision History	52
Appendices	53
Appendix A: Supplies	53
Appendix B: Laboratory Submission Forms	54
Appendix C: Major Foodborne Diseases: General Features	56
Appendix D: Water Testing Fact Sheets	67
Appendix E: Ontario Public Health Standards and Supporting Legislation	81
Appendix F: Full Revision History	85
References	90
Pavision History	101

Introduction

The Public Health Inspector's Guide to Environmental Microbiology Laboratory Testing (the Guide) is designed to support effective public health practice through the provision of timely information on PHO's Laboratory services and expertise.

As an 'evergreen' document the Guide will be continually edited and updated. Any current changes will be noted in revision histories located at the end of each section of the document. For a full list of changes please see Appendix F. If you print a copy of the Guide please check the Public Health Ontario website to ensure you are using the most recent version.¹

Public Health Ontario appreciates the ongoing engagement and generosity of all who are involved in the delivery of environmental health programs and services in Ontario.

Environmental Matrices Analysed by the PHO Laboratory

- Microbiological analyses of foods
- Environmental swab analyses for pathogen specific requests (outbreaks)
- Drinking water analyses for private citizens who rely on private drinking water systems (Total Coliform and *Escherichia coli* only)
- Drinking water analyses requested by boards of health (BOH) (specific microbiological analysis only)
- Ice made from treated water
- Recreational water facilities and public pools (O. Reg. 565/06) /spas (O. Reg. 428/05)
- Public beach water testing
- Water samples suspected of being contaminated with sewage
- Environmental sampling for *Legionella* species (outbreak/cluster investigations)

Note: <u>Consultation</u> with the microbiologist is required for requests that fall outside of the scope of the PHO Laboratory test menu. Board of health staff are encouraged to contact the <u>PHO Laboratory</u> <u>Customer Service Centre</u> with inquiries.

Each PHO Laboratory site holds a drinking water testing license (Ministry of the Environment, Conservation and Parks) for drinking water tests, which are limited to the scope of testing, and accredited by the Canadian Association for Laboratory Accreditation Inc. (CALA) to ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories.²⁻⁴

The PHO Laboratory uses references and methods as guidance documents for testing from agencies such as <u>Health Canada</u> (HC), <u>American Public Health Association</u> (APHA), the <u>AOAC</u> and the US <u>Food and Drug Administration Bacteriological Analytical Manual</u> (FDA-BAM).⁵

Investigations and Analyses Not Available at the PHO Laboratory

The PHO Laboratory does not test environmental samples for the following parameters; please see other resources below for testing options:

- Algal toxin testing (Blue-green algae)
- Chemical/radiological analysis
- Foreign objects
- Fungi
- Infestations (e.g., products containing parasites, insects)
- Floatation tank/isopod or foot throne samples
- Virus detection in food or water
- Botulism investigations
- Giardia and Cryptosporidium in water

Other Resources

- <u>Botulism Reference Service</u>, Health Canada⁶
- <u>Canadian Association for Laboratory Accreditation</u> (CALA) for a list of laboratories accredited with CALA in Ontario
- Canadian Food Inspection Agency (CFIA)
- Ministry of Health, Environmental Health Policy and Programs, Health Protection and Surveillance Policy and Programs Branch (MOH-EH) <u>environmental.health@ontario.ca</u>
- Ministry of the Environment, Conservation and Parks (MECP)
- Ontario Ministry of Agriculture, Food, and Rural Affairs (OMAFRA)
- Standards Council of Canada (SCC) for a list of laboratories accredited with SCC in Ontario
- University of Guelph Laboratory Services
- Public Health Ontario Laboratory <u>Criteria for Acceptance of Patient Specimens</u> webpage and <u>Test Information Sheets</u> for clinical specimens

Legal Considerations

The microbial analyses and acceptable limits outlined for the various parameters and sample type/source are based on current standards, where available, and/or guidelines based on current scientific evidence. References have been included to assist the user in obtaining further information.

These guidelines are intended to serve as a general statement regarding the environmental microbiological testing services available at PHO Laboratory and are of no legal force or effect. If a local board of health considers submitting environmental samples as part of an investigation with overlapping jurisdictions (e.g., federally or provincially inspected food plant combined with retail food outlet), and there is a discrepancy between federal and/or provincial microbial standards/acceptable limits, the local board of health should consult with its legal counsel.

Legal Samples – Collection Instructions

PHO Laboratory must receive prior notification that a legal sample has been collected and will be submitted.

Legal sample is defined as:

any laboratory sample that has been submitted for testing to the PHO Laboratory with an intact legal seal and for which the chain of custody has been maintained.

It is essential to ensure that legal samples are handled according to proper guidelines and adhere to the following instructions:

- Notify <u>PHO Laboratory Customer Service Centre</u> when legal samples are to be submitted.
- Affix and secure an official PHO Laboratory (see <u>Appendix B</u>) and/or a law enforcement agency legal seal to the legal sample such that the sample cannot be opened without breaking the seal. The official legal seal must be initialed and dated by the individual who has affixed the legal seal. The sealed sample must be placed in an appropriate tamper evident container. Contact the <u>PHO Laboratory Customer Service Centre</u> to obtain official legal seals.
- Once the legal sample is received at the laboratory, a Chain of Custody form (e.g., the
 "Relinquished By" section of the water sample requisition) must be signed by the individual
 delivering the samples, and the legal seal verified as intact by both the individual delivering and
 receiving the package. Refer to <u>Instructions for Official Agencies Submitting Water Samples to</u>
 the Public Health Ontario Laboratory for additional details on legal submissions.⁷

Chain of custody is defined as:

an unbroken, documented, chain of events from the time of sample collection, through sample transport, testing and reporting.

Maintaining a proper chain of custody from the time of collection of the sample through transport, testing and reporting ensures the integrity of the legal sample remains intact and supports the validity of the laboratory results in a legal case. Consult with legal counsel and the applicable legislative and regulatory requirements as to when legal seals are required.

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change	
May 6, 2019	 Updated the number of boards of health within the province. Environmental swabs submitted for RRO 1990, Reg. 562: Food Premises, Section 80 will no longer be accepted for testing. All applicable areas have been updated to reflect this change. Updated the Ministry of the Environment, Conservation and Parks (MECP) name and link. 	

Food Sample Analysis

Food - Microbiological Analyses Test Menu

PHO Laboratory performs testing of food samples submitted by Ontario's boards of health (BOH) for microbiological analyses including indicator and pathogenic organisms as well as physicochemical parameters. Food samples may be collected to assist with an outbreak investigation (part of meal, control or follow-up (after remedial action)) or for microbial risk assessment in support of a Hazard Analysis Critical Control Point (HACCP) assessment or routine surveillance. Food samples may also be submitted as a 'complaint' submission (food with off-odour, off-colour, texture, etc.) for laboratory evaluations that are not associated with an investigation of foodborne illness.

Analyses available at the PHO Laboratory for food testing are listed below. Test information including testing frequency, turnaround times and reporting limits, is dependent on the specific analysis. [Click on individual testing links for <u>test directory</u> details – see below.]⁸

INDICATOR BACTERIA

Aerobic Colony Count⁹

Total Coliform¹⁰

Escherichia coli¹¹

Total Gram Negative¹²

Yeast and Mould¹³

FOODBORNE PATHOGENS

<u>Bacillus cereus</u>14

Campylobacter species¹⁵

Clostridium perfringens¹⁶

Escherichia coli O157:H7¹⁷

Listeria monocytogenes18

Salmonella species¹⁹

Shigella species²⁰

Staphylococcus aureus²¹

Staphylococcus aureus toxin²²

Vibrio species²³

Yersinia enterocololitica²⁴

PHYSIOCHEMICAL

pH²⁵

Phosphatase²⁶

Water activity (aw)27

Routine testing of cooked, ready-to-eat foods includes: Aerobic Colony Count, *Bacillus cereus*, Total Coliform, *Escherichia coli*, total gram negative, *Salmonella* species and *Staphylococcus aureus*. Other indicator bacteria and foodborne pathogens will be tested depending on the food type, clinical laboratory information and epidemiological details provided on the <u>Food Bacteriology Requisition</u>.

Physicochemical testing will be performed upon request. For a list of etiological agents associated with foodborne illness, including incubation periods, symptoms, modes of transmission and associated foods, refer to Appendix C: Major Foodborne Diseases: General Features.

Analyses NOT available at PHO Laboratory as part of the for food testing program:

- Cryptosporidium
- Cyclospora
- Giardia
- Hepatitis A
- Norovirus

- Prion disease
- Toxoplasma
- Trichinella
- Ciguatera poisoning
- Paralytic shellfish poisoning

The above listed test requests will be referred to <u>CFIA</u>, Ontario-Area Recall Coordinator (ON-ARC), as per <u>Ontario Foodborne Illness Outbreak Response Protocol (ON-FIORP) - Annex 2</u>. ²⁸⁻²⁹ Additional information such as commercial product details, reason for test request, clinical implications and priority level will be required by CFIA for the test referral. Contact the environmental laboratory through <u>PHO Laboratory</u> <u>Customer Service Centre</u> at 416-235-6556 or toll free 1-877-604-4567 to discuss any other test request not listed above.

Botulism test requests - contact <u>Health Canada's Botulism Reference Service</u> directly at (613) 957-0902 or after-hours at (613) 296-1139.⁶

Food Surveys or Studies

If a board of health would like to initiate a special food survey or study that involves laboratory testing, contact the environmental laboratory program microbiologist through PHO Laboratory Customer
Service Centre at 416-235-6556 or toll free 1-877-604-4567 to discuss the study outline and impact on the laboratory.

Sample Selection

Instructions

- Selection of food items for collection and submission related to a foodborne outbreak should be
 determined based on all available environmental, epidemiological and clinical laboratory
 information related to the outbreak. Information about the clinical pathogen, if known or
 suspected, should be included to facilitate prioritization of food items both collected and tested.
- If unsafe food-handling practices are identified during the inspection of a food premise, selection of food items for a microbial risk assessment is based on the health inspector's discretion. The submission of foods is not a requirement in these circumstances but may provide an additional tool available to support the assessment.

Considerations

- During a foodborne outbreak investigation, the results of an environmental investigation, such
 as a Hazard Analysis Critical Control Point (HACCP) review, and epidemiological analyses, such as
 the calculation of attack rate, will facilitate prioritization of food items that are most likely to be
 involved. If attack rates are unavailable, incubation period, clinical symptoms and epidemiologic
 information about related cases should be used.
- If risk stratification of food sources is not performed by the submitting public health inspector, and there is an abundance of samples submitted for testing, prioritization will be determined by PHO Laboratory Environmental Microbiology section. The communication of all available clinical, epidemiological and environmental information to PHO Laboratory is of great importance in guiding these decisions.
- Avoid submitting raw foods which require processing or cooking before consumption, unless
 they are linked to a specific outbreak or a confirmed etiological agent. If testing of raw foods is
 required, please contact the <u>PHO Laboratory Customer Service Centre</u> at 416-235-6556 or toll
 free 1-877- 604-4567.

Containers for Sampling

Instructions

- Sterile plastic bags with wire closure are the preferred sample container for all food types.
- Food samples found in the original container can be submitted in the original packaging providing that the container is air-tight and will not leak during transit.
- If the original container has been opened or if the packaging is fragile and may be damaged during transit, samples should be transferred from the original container to a sterile plastic sampling bag with wire closure.

Considerations

- If a commercial packaged product has been submitted, record the following product information on the Food Bacteriology Requisition: brand name, lot number, best before date, expiry date, use by date, Universal Product Code (UPC), name of manufacturer/importer/legal agent, unit size, non-intact (open) and intact (closed). Photographs may also be taken and maintained by the public health inspector to capture this information at the time of sample collection which may facilitate a food recall if required. Photographs should not be submitted to the laboratory.
- For foodborne outbreak investigations: when collecting an opened commercial packaged product, collection of a closed sample with the same lot may be useful to determine if the product was contaminated prior to consumer handling. PHO Laboratory will only analyze closed samples in conjunction with a foodborne outbreak investigation and only following consultation and approval by the laboratory. Please consult with the CFIA prior to submission as appropriate.²⁸
- Styrofoam and box-board packaging used for fast-food products are examples of containers that are **not suitable** for microbiological testing. They are not air tight and may leak, potentially resulting in cross contamination of sample(s). These containers **should not** be put inside a sterile sample bag since the outside surfaces of the containers are a potential source of contamination and the sample will be compromised.

Sample Preparation

Instructions

The following steps are recommended to organize and minimize food sample collection time:

- Investigate before sampling to determine a sampling plan.
- Ensure all materials required are available prior to collection.

- Label sample bags using a permanent marker.
- Complete requisition(s) before or after, never during, sample collection.
- After collection, double check to ensure the food sample unique identifier on the sample bag and requisition are legible and match.

Sample Collection

Instructions

- Use aseptic technique at all stages of sample collection. If sterile gloves or tongs are not available, an inverted sampling bag may be used as a glove to collect the sample.
- Collect samples that are representative of each item of food that is being investigated.
- Meals that contain multiple discreet items (e.g. rice and chicken) may be collected individually
 or combined as a single sample. If collected as a single sample, components should be collected
 in approximately equal amounts and will be tested as one sample.
- Foods that contain multiple items (e.g., sandwiches) may be collected as a single sample.
- If an individual component is of greater interest than other components of a meal or mixed food, sample the items individually so they can be tested as individual samples.
- Once collection is completed, close the mouth of the sterile plastic sampling bag with the wire closure, and roll down several times twisting the ends together tightly to prevent leaking during transport.
- Submit a minimum of 200 grams of each sample if possible, as multiple tests may be required. Samples that are less than 10 grams will not be tested.

Considerations

- Failure to follow these instructions will compromise the sample integrity and the laboratory results.
- Separation of meals or mixed foods must be performed at the sampling stage. The laboratory will not separate components of a sample once received in the laboratory.
- Leaking samples will not be accepted by the laboratory. Refer to the Food Laboratory
 Acceptance Criteria section for a complete list of acceptance criteria.

Sample Collection – Large Lots of Food

Most pre-packaged foods are produced in very large lots (a batch or production unit which may be identified by the same code; i.e., lot number, production code or other unique identifiers).

Instructions

• Where possible, mix the same lot of food before sampling.

Considerations

- If a lot of food is larger than 2 kilograms, take five sub-samples of 200 grams each and submit as individual samples for analysis.
- Ensure sub-samples are clearly and distinctly labeled with a unique identifier for tracking purposes on both the sampling bag and requisition.

Sample Collection – Small Quantities of Food (<200 grams)

During some investigations, it may be determined that there is very little sample available for analysis.

Instructions

- In foodborne outbreak investigations, small quantity samples can be processed by the laboratory, but complete analysis may not be possible. Tests performed will be prioritized based on the type of food/suspected etiological agent.
- Samples that are less than 10 grams will not be tested.

Considerations

• If a particular etiological agent(s) is suspected or confirmed, indicate the information on the requisition form and the analysis for this agent(s) will be given priority. In the absence of such information, the laboratory will aim to perform those analyses most likely to provide information useful to the investigation.

Completion of the Food Bacteriology Requisition

The information below outlines the instructions for completing the <u>Food Bacteriology Requisition</u>. Complete instructions are also located on the reverse or page 2 of the requisition and include, submitter information, sample collection, reason for test request, sample details and shipping instructions.

Instructions

• Complete all areas of the <u>Food Bacteriology Requisition</u>. The information included on the requisition will also be provided on the final laboratory report.

Shipping Instructions

- Ship dry foods and other shelf stable foods in a closed container at ambient temperature.
- Ship frozen foods in an insulated container with sufficient ice packs to maintain the frozen state.
- Ship perishable foods in an insulated container with sufficient cold packs to maintain a temperature as close to 4°C as possible. If ice is used, contain the ice in a manner that does not allow water contact with the samples.
- Submit all food samples to the closest PHO Laboratory as soon as possible after collection.

Considerations

Foods that exhibit temperature abuse during shipping may be rejected and not processed. If the
product appears to have undergone a change of state such as from a frozen state to a liquid
state, the product will be rejected.

Food – Laboratory Acceptance Criteria

The accuracy of the test results may be affected by improper collection, handling and/or shipping. Food samples that do not meet the acceptance criteria may be rejected by the laboratory and a new sample would be required to be re-submitted with a newly completed requisition. The acceptance criteria are outlined below.

Description	Acceptance Reason		
Submitter	Sample must be received from an authorized board of health.		
Sample Type	 Food samples are collected to assist with an outbreak investigation (part of meal, control or follow-up (after remedial action)), complaint or for a routine monitoring microbial risk assessment (HACCP or routine surveillance). PHO Laboratory may analyze raw food samples (e.g., raw chicken), which require processing or cooking before consumption, to support a foodborne outbreak investigation with a confirmed etiological agent. 		
Sample Requisition	 The requisition must be completed when it is received at the laboratory including date, time and location of collection. The 		

Description	Acceptance Reason
	 requisition must have a unique identifier that matches the identifier on the sample collection bag. If a sample is received without a requisition it will not be processed; however, the sample will be held for two weeks.
Sample Collection	 The sampling container must not compromise the integrity of the sample. Refer to the Containers for sampling section. The full panel of tests may not be performed if a small amount of food (<200 grams) is received. At least 25 grams of sample is required for each analysis. Please contact the PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567, if collecting a small amount of food.
Sample Transport	 The sample container must be secure to avoid leaks during transport. The sample temperature must meet the requirements listed below when received at the laboratory: Frozen foods should be received frozen (≤ 0°C); if frozen foods have a temperature ≥ 15°C when received at the laboratory, they will not be tested. Foods normally held at refrigeration temperatures should be between 2.0 to 8.0°C when received at the laboratory and must not be at or < 0°C or ≥ 25 °C. Shelf stable foods should be shipped at ambient temperatures and must be between 0°C and 30°C when received at the laboratory.

Food – Reporting and Interpretation of Environmental Microbiological Test Results

Satisfactory microbial limits for a food sample will be dependent on the food preparation, conditions (raw, ready-to-eat, and frozen, refrigerated, shelf stable) and commodity type (e.g., dairy, meat, etc.) in relation to a specific analysis. Refer to Health Canada's Health Products and Food Branch (HPFB).

Standards and Guidelines for Microbiological Safety of Food - <u>An Interpretive Summary</u> for food commodity type details.³⁰

The tables below describe the reporting limit and unsatisfactory level for ready-to- eat foods for indicator bacteria and foodborne pathogens according to Health Canada, and physicochemical testing available at PHO Laboratory. If a board of health would like further interpretation on food results, contact PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567, and provide the sample identification number.

INDICATOR BACTERIA

Testing Type	Reporting Limit (CFU/g) ^a	Unsatisfactory Level (CFU/g) a, c
Aerobic Colony Count	< 1,000 to > 200,000	Category 1: $\geq 10^5$, Category 2: $\geq 10^7$, Category 3: N/A
Total Coliform	< 3 to > 1,100 b	≥ 10 ³
Escherichia coli	< 3 to > 1,100 ^b	≥ 100
Total gram negative	< 1,000 to > 200,000	≥ 10 ^{4 d}
Yeast and mould	< 10 to > 200,000	No reportable limit available

a - CFU (colony forming units)

Unsatisfactory levels of indicator bacteria indicate improper food preparation, handling, storage or holding conditions, temperature abuse, and/or sanitation. Refer to specific test details below for more information.

<u>Aerobic Colony Count</u> (ACC): It is reasonable to expect properly prepared and handled ready-to-eat processed or cooked foods, where no additional food handling was required after food preparation, to have an ACC of $< 10^4$ CFU per gram. An ACC of $> 10^5$ CFU per gram is usually evidence of post processing temperature abuse, inadequate cooking or improper storage conditions. ACC cannot be applied to raw, uncooked, unprocessed foods (e.g., fruits or vegetables), cultured or fermented products (e.g., cheese, yogurt or salami, etc.) since these foods are known and expected to have naturally occurring microbial flora. There are three categories of ACC listed in the table above, which are based on food type and the processing/handling the food has undergone.

Category 1 – cooked foods that do not require handling or processing; i.e., soups, bread, quiche, cooked meat, fish & seafood, vegetables, etc.

Category 2 – cooked foods that require further handling prior to or during the preparation of the final product; i.e., hot dogs, sandwiches, burgers, etc.

Category 3 – foods that have a high ACC due to the normal microbial flora associated; i.e., pitas, potato or pasta salad, salad rolls, fresh fruits, raw vegetables, cheese, yogurt, deli meats, etc.

<u>Total Coliform</u>: Presence of coliforms in a processed, cooked, ready-to-eat food is indicative of inadequate processing or post processing contamination resulting from poor handling of the product. A

b - MPN (most probable number) is used for this test method

c – Microbial guidelines as per Health Canada

d - Microbial guidelines as per Alberta Health Services³²

satisfactory limit cannot be applied to uncooked or raw foods such as salads and salad ingredients, raw fruits and vegetables, raw meats as well as some fermented foods, etc., since coliforms are regularly found in and/or on these foods.³³

<u>Escherichia coli</u>: The presence of <u>Escherichia coli</u> in food is indicative of direct or indirect contamination of the food product with fecal material, thus, enteric pathogens may also be present. The presence may be indicative of inadequate processing, post processing contamination and poor sanitation.

<u>Total gram negative</u>: Total gram negative (TGN) count has the same significance as presence of coliforms. Elevated TGN count should not be detected in cooked or processed ready-to-eat foods. Absence of gram negative organisms in a food sample is another indicator of food safety.

<u>Yeast and mould:</u> The presence of yeast and moulds can cause various degrees of deterioration and decomposition of food, resulting in abnormal odours, flavours and colours. Both yeast and mould can also be added intentionally to a product for a desirable flavour such as mould for blue and brie cheeses and yeast for beer, ciders and wines.

FOODBORNE PATHOGENS

Testing Type	Reporting Limit (CFU/g) ^a	Unsatisfactory Level (CFU/g) a, b
Bacillus cereus	< 100 to > 200,000	≥ 10 ⁴ potentially hazardous
<u>Campylobacter</u> species	Not Detected / Detected	Detected
Clostridium perfringens	< 100 to > 200,000	≥ 10 ⁴ potentially hazardous
Escherichia coli O157:H7	Not Detected / Detected	Detected
Listeria monocytogenes	Not Detected / Detected	Detected
Salmonella species	Not Detected / Detected	Detected
<u>Shigella species</u>	Not Detected / Detected	Detected
Staphylococcus aureus	< 100 to > 200,000	≥ 10 ⁴ potentially hazardous
Staphylococcus aureus toxin	Not Detected / Detected	Detected
<u>Vibrio species</u>	Not Detected / Detected	Detected
Yersinia enterocolitica	Not Detected / Detected	Detected

a - CFU (colony forming units)

b – Microbial guidelines as per Health Canada³⁰

Foodborne pathogens should not be detected from a ready-to-eat food and would indicate improper food preparation, poor handling, inadequate sanitation practices or possible contamination from a positive food source or an ill food handler. <u>Appendix C</u>: Major Foodborne Diseases: General Features provide additional information regarding the organisms incubation period, clinical symptoms, mode of transmission and associated foods.

PHYSICOCHEMICAL

Testing Type	Reporting Limit
рН	1.00 to 14.00
<u>Phosphatase</u>	Not Detected / Detected
Water activity (a _w)	0.000 to 1.000 IU/L

<u>pH:</u> The pH is a function of the hydrogen ion concentration in the food and is a measure of food acidity and alkalinity, which varies with food type. pH is measured on a scale of 1 to 14, with 7 being neutral, lower numbers indicating greater acidity and higher numbers indicating greater alkalinity. pH is just one factor that influences bacterial growth, and can interact with the following to inhibit growth of pathogens and other organisms: water activity (a_w), salt, redox potential, preservatives and temperature. As the pH decreases, a lower temperature is needed to inactivate microbes. The phodes with an equilibrium pH \leq 3.7 will not support the growth of bacterial foodborne pathogens, regardless of the a_w. However, if the pH of the food (jarred/canned products) is \leq 4.6, microorganisms are inhibited providing the a_w is \leq 0.85. According to Public Health Ontario's Home Canning Literature Review, June 2014, a hermetically sealed acidic canned or jarred food with an equilibrium pH \leq 4.6 that has gone through sufficient heat treatment to eliminate vegetative microorganisms (bacteria, yeasts and molds), will not support the germination and growth of botulinum spores and production of the toxin.

<u>Phosphatase</u>: The enzyme alkaline phosphatase (ALP) is associated with the cream or fat globules of raw milks. It is inactivated below the detection limit of conventional methods following complete pasteurization. The detection of residual alkaline phosphatase indicates a reduction in pasteurization temperature of at least 1.5° C, a 5 minute reduction in holding time, or the presence of $\geq 0.3\%$ raw milk. Since ALP is more heat stable than most pathogens (which may be present in milk), it is used as an indicator of pasteurization; however, a negative ALP test does not guarantee that the product is pathogen free.³⁶

<u>Water activity (a_w)</u>: The availability of water for microbial, enzymatic or chemical activity determines the shelf life of foods. This water availability is measured as water activity (a_w). Water activity is the ratio of water vapor pressure of the food substrate to the vapor pressure of pure water at the same temperature. Water activity is measured on a scale of 0 to 1; pure water has an a_w of 1.00 and the a_w of a completely dehydrated food is 0.00. For most foods to be considered safe to store at room

temperature, they should have an $a_w \le 0.85$. Bacterial foodborne pathogens are inhibited at this aw; however, for those microorganisms that can cause spoilage (e.g., yeasts and moulds), the a_w should be < 0.60 to prevent food spoilage by inhibiting growth of these organisms.³⁷

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the Public Health Ontario website to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change	
May 6, 2019	• Updated the number of boards of health within the province.	

Environmental Swabs

Swab – Microbiological Analysis Test Menu

PHO Laboratory performs microbiological pathogen specific testing of environmental swabs to support foodborne outbreak investigations. Environmental swabs processed during foodborne outbreak investigations will be tested for the presence or absence of a specified clinically confirmed bacterial pathogen from a physical environment (e.g., cutting board, meat slicer, etc.).

Sampling instructions for environmental swabs are provided below. Test information including testing frequency, turnaround times and reporting limits, is dependent on the specific analysis. [Click on individual testing links for <u>test directory</u> details - <u>Salmonella spp.</u>, <u>Listeria monocytogenes</u> and <u>E. coli</u> O157:H7.]^{8,19,18,17}

Sample Selection

Instruction for sample selection of swabs for an outbreak investigation:

To identify the source of the pathogen, selection of sampling sites for the investigation should be informed by environmental and epidemiological factors. Random sampling may delay the identification of the source of the outbreak. Samples deemed as most likely to be positive should be sampled and submitted first.

PHO Laboratory should be notified prior to sample submission to allow for preparation of materials required for testing. Please contact the laboratory by calling PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567.

Containers for Sampling

Instructions

- Environmental swabs are available from the PHO Laboratory warehouse and can be ordered by calling the <u>PHO Laboratory Customer Service Centre</u> at 416-235-6556 or toll free 1-877-604-4567. The swab vial contains neutralizing solution as transport medium. This medium is able to neutralize quaternary ammonium compounds and phenolic disinfectants.
- Refer to the specific foodborne pathogen test directories, i.e., <u>Salmonella spp.</u>, <u>Listeria</u> monocytogenes or E. coli O157:H7 for additional details. ^{19,18,17}

Each sample and its accompanying requisition must be labeled with a unique identifier.

Sample Preparation

The following steps are recommended to organize and minimize sample collection time:

- Investigate before sampling to determine a sampling plan.
- Ensure all materials required are available prior to collection.
- Label sample container using a permanent marker.
- Complete requisition(s) before or after, never during, sample collection.
- After collection, double check to ensure the sample unique identifier on the sample container and requisition are legible and match.

Sample Collection

Instructions

- Using aseptic technique, open the sterile swab container; grasp the handle end of the swab.
- Remove the swab aseptically being careful not to touch any portion that might be inserted into the vial.
- Moisten the swab head and press out the excess solution against the interior wall of the vial with a rotating motion.
- Hold the swab handle to make a 30 ° angle contact with the surface. Rub the surface slowly and thoroughly over approximately 50 cm² of the surface three times, reversing direction between strokes.
- After the area has been swabbed, using aseptic technique, return the swab to the vial.
- Complete all sections of the <u>Environmental Microbiology Investigation Requisition</u> form with appropriate information to ensure accurate tracking of the sampling and complete all fields as outlined on page 2 of the form.
- Ensure samples are submitted to the laboratory under refrigeration conditions (2 8 °C) in a manner that will not leak during transport.

Considerations

- Leaking swab vials will be rejected by the laboratory. Ensure the vial is closed tightly to prevent leakage resulting in sample rejection.
- Absence of a unique identifier linking sample(s) to a sampling location on a requisition will result in rejection of sample(s) by the laboratory.
- For an investigation where sampling areas are smaller and physically different, a new swab should be used for each different area.
- For an investigation where the sample area is large, multiple sections may be swabbed with one swab. Sampling details should be accurately recorded to ensure traceability when laboratory results are received to ensure appropriate interpretation of results.

Environmental Swab – Laboratory Acceptance Criteria

The accuracy of the test results may be affected by improper collection, handling and/or shipping. Swab samples that do not meet the acceptance criteria may be rejected by the laboratory and a new sample may be re-submitted with a newly completed requisition. The acceptance criteria are outlined below.

Description	Acceptance Reason
Submitter	 Swab samples must be received from an authorized board of health submitter.
Sample Type	 Swabs will only be accepted for testing to support a foodborne outbreak investigation.
Sample Requisition	 The requisition must be completed upon receipt at the laboratory, including date and location of collection. If a sample is received without a requisition it will not be processed; however, the sample will be held for 24 hours. The requisition must have a unique identifier that matches the unique identifier on the swab collection vial.
Sample Collection	 Environmental swabs are available from the PHO Laboratory warehouse and can be ordered by calling PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567. Unapproved swabs will not be accepted.
Sample Transport	 The sample container must be secure to avoid leaks during transport. The sample temperature must meet the requirements listed below upon receipt at the laboratory:

Description	Acceptance Reason	
	Swabs should be shipped with cold packs to maintain a	
	temperature between 2.0 to 8.0 °C during transport to the	
	laboratory.	
	 Swabs received frozen or greater than room temperature 	
	may be rejected and not processed.	

Environmental Swab – Reporting and Interpretation of Environmental Microbiological Test Results

The following tables outline the reporting limits and interpretation of the environmental swab analyses. Information such as collection conditions, time and date, and sanitation details should be considered when interpreting swab results.

FOODBORNE PATHOGENS

Testing Type	Reporting Limit	Unsatisfactory Level
Escherichia coli O157:H7	Not Detected / Detected	Detected
<u>Listeria monocytogenes</u>	Not Detected / Detected	Detected
Salmonella species	Not Detected / Detected	Detected

Foodborne pathogens should not be detected from a ready-to-eat food preparation surface swab and would indicate poor food handling, inadequate work surface sanitation practices or possible contamination from an ill food handler. <u>Appendix C</u>: Major Foodborne Diseases: General Features provide additional information regarding the organisms incubation period, clinical symptoms, mode of transmission and associated foods.

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the Public Health Ontario website to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change
August 8, 2021	 Updated the swab requisition details from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition.

Water Sample Analysis

Water – Microbiological Analysis Test Menu

PHO Laboratory performs microbiological testing of water samples for the detection of bacterial indicators of contamination and specific pathogens in water. Each water sample is tested for specific indicators according to the water source.

If testing is required on a water source that is not listed in the PHO Laboratory water testing menu below, or testing is needed to identify specific etiological agents in water related to a laboratory confirmed clinical case or outbreak investigation, boards of health are asked to <u>consult</u> with the microbiologist, or designate, overseeing the water testing program prior to sample submission.

Sample collection requirements, sample handling, shipping conditions, test information including testing frequency, turnaround times and reporting limits are dependent on the specific water source. Samples will not be processed if the requisition is not completely and accurately filled in when received at the laboratory; and a new sample and completed form will be required to be submitted. [Click on individual testing links for test directory details – see below.]⁸

Drinking water (including bottled water)³⁸

<u>Drinking water – Private citizen</u>³⁹

Ice - Treated⁴⁰

Public Beach Water⁴¹

Recreational Water Facilities, Public pools/Spas⁴²

Suspected Sewage Contamination – Water⁴³

For specific etiological agent requests, please contact the <u>PHO Laboratory Customer Service Centre</u> at 416-235-6556 or toll free 1-877-604-4567.

For general information related to sample collection, labelling, storage, transportation and chain-of-custody of drinking water samples that fall under the Drinking Water Systems Regulation (Ontario Regulation 170/03) and Small Drinking Water Systems (Ontario Regulation 319/08), refer to the Ministry of Environment, Conservation and Park's <u>Practices for the Collection and Handling of Drinking-Water Samples</u>.

Important Information Related to Drinking Water Analysis

Type of Drinking Water System

The type of drinking water system identified on the requisition triggers specific notification requirements. Laboratories have an obligation to report test results as outlined in section 18 of the <u>Safe Drinking Water Act, 2002</u> for drinking water samples submitted under an Ontario drinking water regulation <u>O. Reg. 170/03 Drinking Water Systems</u> or <u>O. Reg. 319/08 Small Drinking Water Systems</u>. 45-47 The details are outlined in the <u>Drinking Water – Official Agencies</u> Test Information Sheet. 38

Public health inspectors may collect drinking water samples in accordance with the following sections of the Health Protection and Promotion Act (<u>HPPA</u>): Section 13, Section 41, Subsections (1), (2), (4), and (5) which include drinking water systems that are not captured by a specific regulation(s), but where the public has access to the water, and are the subject of an investigation under the <u>Health Protection and Promotion Act</u> (HPPA).⁴⁸

These may be considered 'drinking water systems' in which drinking water samples must be analyzed at a licensed laboratory using licensed drinking water testing methods, unless the example is entirely a plumbing system.

The table below provides examples of drinking water systems that fall outside the scope of the <u>Safe</u> <u>Drinking Water Act</u>, <u>2002</u>, <u>Drinking Water Systems Ontario Regulation 170/03</u> and the <u>Health Protection</u> and <u>Promotion Act</u>, <u>Small Drinking Water Systems Regulations O. Reg. 319/08</u>. ⁴⁵⁻⁴⁷ Consult with Ministry of Health, Environmental Health Policy and Programs, Health Protection and Surveillance Policy and Programs Branch (MOH-EH) for direction on drinking water system classifications.

Completing the requisition and submitting the sample according to <u>Instructions For Official Agencies</u> <u>Submitting Water Samples to the Public Health Ontario Laboratory</u>⁷ will help to ensure samples are accepted.

EXAMPLES OF DRINKING WATER THAT FALL OUTSIDE ONTARIO DRINKING WATER REGULATIONS

Situation	Example
Drinking water supplies that	 Seasonal campsite with fewer than 6 water connections,
do not meet the specific	camps in unorganized territory, certain recreational camps,
criteria of the Safe Drinking	rental cabins with fewer than 6 units using the same water
Water Act 2002 or Health	supply, etc. These are camps in which buildings are used to
Protection and Promotion	accommodate five or more employees who are employed in
	mining work, lumbering work or any other labour work in

Situation	Example
Act drinking water regulations	territory without municipal organization (R.R.O. 1990, Regulation 554). 49 Raw water that it is not for human consumption and is not distributed to the public.
Facilities in the following circumstances	 Serviced by private drinking water systems that require health unit approval for licensing or funding approval: Funeral homes, group and boarding homes, homes for special care, migrant farm workers, etc.
Food premises in the following circumstances	 Small food take-out and convenience stores where there is no direct access to drinking water by the patrons, temporary food service facilities, farmers markets, etc. Special event food premises or mobile food vending vehicles
Outbreak investigations	 Testing of a drinking water system as part of an outbreak investigation (e.g., E. coli O157:H7, Legionella, Salmonella, etc.).
Private individual well water supplies	 Enteric disease follow-up Assisting with addressing/identifying water safety concerns Complaints Landlord and tenant matters whereby samples are taken of the private water supply for the issuance of a section 13 order for a violation of section 20 under the HPPA⁴⁸ Treatment devices for samples taken to assess the ability of a newly installed treatment process at a premise
Water Containers	 Water haulage trucks (testing may be required as part of an inspection) Private building water coolers Commercial bottled water/water bottling facilities Vending machines that dispense water, i.e. "Bottle your own water" Cisterns not supplying systems regulated under the <u>Safe</u> <u>Drinking Water Act, 2002</u> or the <u>Health Protection and</u> <u>Promotion Act</u> 45,48

Chain of Custody

The chain of custody must be maintained for drinking water samples collected under Ontario drinking water regulations (i.e., O. Reg. 170/03 Drinking Water Systems or O. Reg. 319/08 Small Drinking Water Systems). 46-47

For drinking water indicated as "HPPA O. Reg 319/08" or "SDWA O. Reg 170/03", there are two ways to maintain the chain of custody:

f the person collecting the sample cannot deliver the ample directly to the laboratory: Place a signed and dated "Regulated Water Seal" (refer to appendix B: Laboratory Submission Forms) up and over the ap of the sample in such a way that the two ends of the eal are attached to the shoulder of the bottle, and the
nidpoint of the seal is stuck to the top of the cap. The sample is sent to the laboratory with a completed equisition according to the sample transport instructions in the test directory — test information sheet.8
riple is delivered directly to the laboratory by the submitter, at the following fields in the "For Regulated Drinking Water Samples" section of the requisition: The printed name of the person relinquishing the sample to aboratory personnel The signature of the person relinquishing the sample Date and time the sample(s) is relinquished

Compliance versus Audit Samples

Owner/operators of O. Reg. 170/03 or O. Reg. 319/08 drinking water systems are required to follow their respective directive as determined by the outcome of the risk assessment performed on the drinking water system and, if sampling is required, use a licensed commercial laboratory for their testing. 46-47

The exception is O. Reg. 319/08 small drinking water systems that supply Ontario Parks; routine compliance samples are analyzed by the PHO Laboratory.⁴⁷

Boards of health may sample regulated drinking water systems as part of an audit; however, they should refrain from routine monitoring of these drinking water systems.

Drinking Water Sample Hold Time

Drinking water must be analyzed in accordance with a test method listed in the Ministry of the Environment, Conservation and Parks (MECP) "Protocol of Accepted Drinking Water Testing Methods" according to the Drinking Water Testing Services regulation O. Reg. 248/03. 50-51

PHO Laboratory uses a modified versions of the MECP methods; namely, E3407 Membrane Filtration Method Using DC Agar for the Simultaneous Detection and Enumeration of Total Coliforms and *Escherichia coli* in Drinking Water, and E3408: The Spread Plate Method for the Enumeration of Aerobic, Heterotrophic Bacteria in Drinking Water. These methods include a sample hold time of 48 hours. ⁴⁶⁻⁴⁷ For results to be suitable for the purposes of the *Safe Drinking Water Act, 2002*, samples must be tested within the hold time. ⁴⁵

Etiological Agent Testing

A Public Health Inspector or designate under the authority of the Medical Officer of Health may collect water samples for identification of an etiological agent (e.g., *Legionella* spp., *Salmonella* spp., etc.) from drinking water systems if required for investigational purposes or outbreak situations. On the <u>Bacteriological Analysis Water - Multiple Sample Requisition for Official Agencies</u>, these samples must be identified as "Non-Regulated", "HPPA Regulated Premises", or "Private Residence", and marked as an "Outbreak Investigation" with the outbreak/investigation number, etiological agent and clinical case identification (e.g., iPHIS, PHO Laboratory barcode identification). Alternatively, complete all sections of the <u>Environmental Microbiology Investigation Requisition</u> with appropriate information to ensure accurate tracking of the sampling and complete all fields as outlined on page 2 of the form.

Drinking Water Submissions under a Provincial Officer's Order or Approval (Directive)

For samples requiring analysis for a bacteriological health related parameter not included in <u>Ontario Drinking Water Quality Standards Ontario Regulation 169/03</u>;⁵² e.g., Heterotrophic Plate Count (HPC) or Background, and submitted under a Provincial Officer's Order or Approval (Directive),⁵³ from a regulated drinking water system, please contact the <u>PHO Laboratory Customer Service Centre</u> and ask to speak to the Program Coordinator Drinking Water Testing or designate.

Drinking Water Research

The <u>Drinking Water Testing Services regulation O. Reg. 248/03</u> outlines the requirements related to drinking water quality research (e.g., scientific studies, including special audits or surveys to assess drinking water quality).⁵¹

If board of health personnel wish to engage in a study involving drinking water (e.g., to identify trends associated with ground water contamination in various geographic areas, etc.) contact PHO Laboratory Customer Service Centre and ask to speak to the Program Coordinator, Drinking Water Testing.

Research projects must be reviewed by the microbiologist overseeing the PHO Laboratory Environmental Microbiology program and approval must be granted by the PHO Laboratory Research and Project Review Committee before research begins.

Water – Laboratory Acceptance Criteria

The accuracy of the test results may be affected by collection, handling and shipping. Submit water samples in a timely manner to the PHO Laboratory. Water samples that do not meet the acceptance criteria may not be accepted by the laboratory. A new sample may be submitted with a new requisition. The acceptance criteria are outlined below.

Description	Acceptance Reason
	 Official Agencies: Local board of health or an Ontario Park under the Ministry of Natural Resources and Forestry (MNRF). Requests for testing by any other organization (e.g., other public health partners or Ontario Ministries) will require
Submitter	 approval from the Director Quality, Regulatory Affairs and Customer Service and microbiologist overseeing the water testing program, or designate, prior to submission. Private drinking water systems (e.g., well water): Private citizen submissions of samples collected from a "single household" and serves fewer than 6 private residences but not the public.
	 The water sample must be collected from within the province of Ontario. The following types of water sources are analyzed at the PHO Laboratory. For specific analyses, refer to the PHO Laboratory test directory – test information sheet:⁸
Sample Type	Local board(s) of health:

Description	Acceptance Reason
	 MNRF Ontario Parks: Drinking Water – Official Agency Public Beach Water Private Citizen Drinking water – Private citizen Refer to the acceptance criteria
Sample Requisition	 The sample must be received with a requisition that is complete and accurate. A unique identifier is used to link the sample with the requisition: For Official Agency water samples, the sample barcode or barcode number must be present on both the requisition and the sample when it is received at the laboratory. For treated ice samples, a unique identifier is written on the sample bag and the requisition; they must match. For Private Citizen drinking water samples, a unique identifier, e.g., the submitter's name can be used as the unique identifier to link the requisition and the sample; however, samples are accepted without a unique identifier when the sample is received at the laboratory with the requisition secured around the bottle. Water samples for identification of an etiological agent (e.g., Legionella spp., Salmonella spp., etc.): These samples must be identified as "Non-Regulated", "HPPA Regulated Premises", or "Private Residence", and marked as an "Outbreak Investigation" with the outbreak/investigation number, etiological agent and clinical case identification (e.g., iPHIS, PHO Laboratory barcode identification).
Sample Collection	 Water samples must be collected in PHO Laboratory water collection bottles. Drinking water samples must be collected from the cold water tap and not the hot water tap. Treated ice samples must be collected in a sterile plastic bag with wire closure.

Description	Acceptance Reason
	 Sufficient amount of sample for the analyses is required (i.e., 200 mL unless otherwise specified in the test directory – test information sheet).⁸ Drinking water samples must be aesthetically acceptable; i.e., clear, colourless and odourless, and does not contain debris or substances that interfere with the test. For Private Citizen drinking water submissions, only requisitions with a version number of 3743-44 (01/01) or newer are accepted.
Sample Transport	 The sample container must be secure and not broken or damaged when received at the laboratory. Samples should be stored and transported at temperatures between 2 and 8 °C following collection. Drinking water samples must be less than 25 °C but not frozen when received at the laboratory. For details on sample transport for public beach, recreational water facility or public pool/spa refer to the test directory – test information sheet.8
Sample Hold Time	 Drinking water samples must be sent to the laboratory as soon after collection as possible since they must be tested within the hold time for the test as specified in the test directory – test information sheet for the analysis required. Non potable water (e.g., public beach, public pool/spa, recreational water facility, sewage, etc.) must be tested within one calendar day of collection.

Water – Reporting and Acceptable Limits of Environmental Microbiological Test Results

The following tables outline the reporting and acceptable limits for potable and non-potable water testing. Refer to <u>Appendix D</u>: Water Testing Fact Sheets for interpretation and additional details on Background or Non-target bacteria, *Campylobacter* spp., Coliforms, *Escherichia coli*, Heterotrophic Plate count (HPC), *Legionella* spp., *Pseudomonas aeruginosa*, *Salmonella* spp. and *Staphylococcus aureus*.

POTABLE WATER

OFFICIAL AGENCY – POTABLE / ICE-TREATED

Testing Type	Reporting Limit	Acceptable Limit
Background or Non-target bacteria (CFU / 100 mL)	0 to *NDOGN/ NDOGT	At the discretion of the MOH. Only reported on Regulated drinking water samples submitted under a "Provincial Officer's Order" or "Approval (Directive)".
Total Coliform ^a (CFU / 100 mL)	0 to > 80	0
Escherichia coli ^a (CFU / 100 mL)	0 to > 80	0
Heterotrophic Plate Count (HPC) (CFU / 1 mL)	< 10 to > 3.0 x 10 ³ , **NDOGHPC	< 500

^{*}NDOGN - No Data Overgrown with non-target/NDOGT - No Data Overgrown with target

OFFICIAL AGENCY – BOTTLED WATER

Testing Type	Reporting Limit	Acceptable Limit
Total Coliform ^a (CFU / 100 mL)	0 to > 80	0
Escherichia coli ^a (CFU / 100 mL)	0 to > 80	0
Heterotrophic Plate Count (HPC) (CFU / 1 mL)	< 10 to > 3.0 x 10 ³ , *NDOGHPC	< 500
Pseudomonas aeruginosa (CFU / 100 mL) (CFU / 50 mL) ^b	0 to > 100 < 2 to > 200 b	0
Staphylococcus aureus (CFU / 100 mL) (CFU / 50 mL) ^b	0 to > 100 < 2 to > 200 b	0

^{*}NDOGHPC - No Data Overgrown heterotrophic plate count

^{**}NDOGHPC - No Data Overgrown heterotrophic plate count

PRIVATE CITIZEN – PRIVATE DRINKING WATER SYSTEMS

Testing Type	Reporting Limit	Acceptable Limit
Total Coliform ^a (CFU / 100 mL)	0 to > 80	≤ 5
Escherichia coli ^a (CFU / 100 mL)	0 to > 80	0

^a An overgrown sample; that is, crowding and/or confluent and/or non-identifiable microbial growth on a plate, is considered an adverse result. When there is an overgrown condition and there is no evidence of Total Coliform or *E. coli*, NDOGN - No Data Overgrown with non-target is reported. The water may be unsafe to drink. When there is an overgrown condition and there is evidence of Total Coliform and/or *E. coli*, NDOGT - No Data Overgrown with target is reported. The water is unsafe to drink.

NON POTABLE WATER

OFFICIAL AGENCY - PUBLIC BEACH WATER

Testing Type	Reporting Limit	Acceptable Limit	
Escherichia coli	< 10 to > 1.0 x 10 ³	≤ 200	
(CFU / 100 mL)	\ 10 to / 1.0 \ 10	2 200	

OFFICIAL AGENCY – SUSPECTED SEWAGE CONTAMINATION

Testing Type	Reporting Limit	Acceptable Limit
Escherichia coli	< 1.0 x 10 ³ to	Not Applicable
(CFU / 100 mL)	> 1.0 x 10 ⁶	ног Аррисавіе

OFFICIAL AGENCY - RECREATIONAL WATER FACILITIES, PUBLIC POOLS/SPAS

Testing Type	Reporting Limit	Acceptable Limit
Background or Non-target bacteria (CFU / 100 mL)	Only reported if NDOGN/NDOGT	Not Applicable
Total Coliform ^a (CFU / 100 mL)	< 2 to > 160	< 2

^b When 50 mL is analyzed.

Testing Type	Reporting Limit	Acceptable Limit
Escherichia coli ^a (CFU / 100 mL)	< 2 to > 160	< 2
Heterotrophic Plate Count (HPC) (CFU / 1 mL)	< 10 to > 3.0 x 10 ³ , *NDOGHPC	Refer to Appendix D: Water Testing Fact Sheets
Pseudomonas aeruginosa (CFU / 100 mL) (CFU / 50 mL) ^b	0 to > 100 < 2 to > 200 b	< 10
Staphylococcus aureus (CFU / 100 mL) (CFU / 50 mL) ^b	0 to > 100 < 2 to > 200 b	< 50

^a An overgrown sample; that is, crowding and/or confluent and/or non-identifiable microbial growth on a plate, is considered an adverse result. When there is an overgrown condition and there is no evidence of Total Coliform or E. coli, NDOGN - No Data Overgrown with non-target is reported. When there is an overgrown condition and there is evidence of Total Coliform and/or E. coli, NDOGT - No Data Overgrown with target is reported.

^b When 50 mL is analyzed.

^{*}NDOGHPC - No Data Overgrown heterotrophic plate count

Water Status Codes and Associated Interpretation Statements for Drinking Water Reports

Status code	Interpretation
0	Not yet tested
1	No significant evidence of bacterial contamination
2	Significant evidence of bacterial contamination
3	UNSAFE TO DRINK, evidence of fecal contamination
4	No Data: Overgrown with Non- Target
5	The sample was taken from an unprotected source, such as a lake or river
7	The sample was too old
8	The bottle was received broken or damaged
10	Aesthetically unacceptable - The appearance or odour makes the sample unacceptable as drinking water
11	The requisition was received separated from the sample bottle
12	The sample was not collected in the proper bottle
13	The sample was received frozen

Reports		
Status code	Interpretation	
14	The sample was collected from a hot water tap	
15	Insufficient sample was submitted	
17	The sample leaked in transit	
18	The sample was collected from a source located outside of Ontario	
19	The sample was collected from a municipal water supply	
20	Insufficient information was supplied on the sample	
21	Not tested for some other reason	
22	No Data: Overgrown with Target	
24	Sample received was very warm upon receipt	
25	Bottled water submitted	
26	Interfering substances in the sample	
27	Unique identifier missing	
28	Outdated collection kit received	
29	Unauthorized submitter	
30	Laboratory error	
50	Technical difficulties with the barcode	

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change	
August 8, 2021	 Added the option to submit water samples using <u>Environmental</u> <u>Microbiology Investigation Requisition</u> for etiological agent testing. 	

Legionella Investigations

Legionella - Microbiological Analysis Test Menu

PHO Laboratory performs environmental testing for *Legionella* for investigations of confirmed culture clinical cases or outbreaks of legionellosis. Clinical case information must be provided to the laboratory when environmental testing is requested and prior to sample collection. Sample collection requirements, sample handling, shipping conditions, test information including testing frequency, turnaround times and reporting limits are dependent on the specific matrix. [Click on individual testing links for <u>test directory</u> details - *Legionella* (water) and *Legionella* (swab).]^{8,54-55}

The <u>WHO</u> outlines Legionella infection by category; community acquired, travel associated and nosocomial infection in relation to environmental sampling sites:⁵⁶

Category	Community Acquired	Travel Associated	Nosocomial
Modes of transmission	Inhalation of contaminated aerosol	Inhalation of contaminated aerosol	Inhalation of contaminated aerosol, aspiration, wound infection
Sources of Legionella	Cooling towers; hot and cold-water systems; spa pools, thermal pools, springs; humidifiers; domestic plumbing; potting mixes and compost	Cooling towers; hot and cold-water systems; spa pools, thermal springs and pools; humidifiers	Cooling towers; hot and cold-water systems; spa pools, natural pools, thermal springs; respiratory therapy equipment; medical treatment
Reservoir of Legionella	Industrial sites, shopping centres, restaurants, clubs, leisure centres, sports clubs, private residences	Hotels, cruise ships, camp sites, shopping centres, restaurants, clubs, leisure centres, sports clubs	Hospitals, medical equipment

Sample Selection

Instructions

- To identify the source of the Legionella, selection of sampling sites for the investigation should
 consider both environmental and epidemiological factors. The CDC provides an example of the
 Legionella Environmental Assessment Form which can be used to document details on the water
 systems, and information on developing a sampling plan.⁵⁷
- Legionella outbreaks occur because of aerosolization and inhalation of high numbers of
 Legionella. Investigators should try to determine where in a water system there is stagnation of
 water and subsequent aerosolization of stagnant water; for example any dead ends in a water
 system may allow Legionella to grow to high numbers. L. pneumophila survived and multiplied
 in water at temperatures between 20 °C and 50 °C, with an optimal temperature range of 32–42
 °C
- Samples from locations based on the environmental assessment and most likely to be positive and associated with the clinical case exposure history should be collected and submitted.

The <u>WHO</u> also provided a list of environmental sampling sites which was adapted from <u>Barbaree JM</u>, <u>et al. 1987</u>. ^{56,58} This information can be found on the table below and it can be used as a guide for selection of sampling sites, but it is not a list of mandatory sampling sites.

Environmental Sites	Examples
Potable water outside or on the boundary of the health- care facility property	Treatment plant (raw and refined water), guardhouse or other facility if water is not fed from health-care facility, fire hydrants
General potable water system for healthcare facility	Incoming water pipe(s), water softener (pre and post), preheater (discharge side), primary heater (discharge side), circulating pumps, holding tanks (cold water, discharge side), expansion tank for hot water, back drain on sprinkler system(s), fire line where it branches off main system, water used for respiratory therapy equipment, outlets in patients' rooms
Air compressor system	Vacuum water source Positive pressure equipment side: Condensate from tank(s), water separator(s) directly off compressors, water source(s) near air intake(s), air samples where patients were ill with legionellosis

Environmental Sites	Examples
Potable water final distribution outlets	Haemodialysis water source, before or after demineralizer
Intensive care units	Respiratory therapy (patient's room), cardiac, services with different geographical locations, ice-maker (entry water) and ice
Air-conditioning system	Air handling unit serving area where disease occurred Cooling towers: return from heat exchanger to water (spray/trough and gutter) distribution or pond (sump), water supply
Hot tubs	Pool and balance tank (if fitted), jets and pipes (swab)
Other	Decorative fountain, creeks, ponds, sites of stagnant water

Considerations

- Ideally, sampling should be performed before disinfection of the water system.⁵⁷
- Dead-ends in the water system may be disconnected where pipes are capped instead of completely removed. Any fixtures and water lines that are not used on a regular basis may allow *Legionella* to grow to high numbers.
- Random sampling may delay the identification of the source of the outbreak. The most obvious sampling sites would be water sources in close proximity to the suspected clinical case exposure.

The table below describes types of environmental sampling sites with typical outbreak incubation period pattern and implicated *Legionella* species:⁵⁶

Sampling Site	Disease Outbreaks	Commonly Implicated Legionella species
Cooling water systems	Rapid onset over wide area, resolve within incubation period	Predominantly <i>L. pneumophila</i> sg 1
Hot and cold-water systems	Low numbers of cases over prolonged periods	L. pneumophila sg 1, 2, 4, 6, 12, L. micdadei, L. bozemanii, L. feeleii and others
Hot tubs, natural spa pools, thermal springs	Rapid onset confined to users and those in close proximity	L. pneumophila sg 1, L. micdadei,

Sampling Site	Disease Outbreaks	Commonly Implicated Legionella species
		L. gormanii, L. anisa
Humidifiers, respiratory equipment	Low numbers over prolonged periods. Rapid onset confined to users and those in close proximity	L. pneumophila sg 1, 3, and others
Potting mixes, compost	Low numbers of cases over prolonged periods	Exclusively <i>L. longbeachae</i>

Containers for Sampling

Instructions

- PHO Laboratory water collection bottles and environmental swabs are available from the PHO Laboratory warehouse and can be ordered by calling the PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567. PHO Laboratory water collection bottles contain sodium thiosulfate which is a dechlorination agent used to remove any residual halogen such as chlorine and prevents continuation of bacterial action during sample transit. The swab vial contains neutralizing solution as transport medium. This medium can neutralize quaternary ammonium compounds and phenolic disinfectants. [Click on individual testing links for test directory details for information on obtaining supplies from PHO Laboratory for Legionella outbreak investigations Legionella (water) Legionella (swab)].8,54,55
- Prior to collecting environmental samples from sources other than water or environmental swabs, contact the environmental laboratory through PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877-604-4567.
- Detailed instructions are also available in the reference document <u>CDC Sampling Procedure and</u>
 Potential Sampling Sites.⁵⁹

Sample Preparation

The following steps are recommended to organize and minimize sample collection time:

• Investigate before sampling to determine a sampling plan. Swab samples should be collected first to capture any biofilm from the shower head or aerator. Then collect a water sample slowly into the PHO Laboratory water bottle to collect any dislodged material.

- Ensure all materials, including appropriate personal protective equipment, are available prior to collection.
- Label sample containers using a permanent marker.
- Complete requisition(s) before or after, never during, sample collection.
- After collection, double check to ensure the sample unique identifier on the sample container and requisition are legible and match.

Sample Collection - Swabs

Instructions

- Use one swab per site, e.g., jets or shower head.
- Refer to the test directory link for sample collection details Legionella (swab).55
- Complete all sections of the <u>Environmental Microbiology Investigation Requisition</u> form with appropriate information to ensure accurate tracking of sampling and complete all fields as outlined on page 2 of the form.

Considerations

- Leaking swab vials will be rejected by the laboratory. Ensure the vial is closed tight to prevent leakage.
- Absence of a unique identifier linking sample(s) to a sampling location on a requisition will result in rejection of sample(s) by the laboratory.

Sample Collection - Water

Instructions

- Using aseptic technique, collect the water sample in PHO Laboratory water collection bottles.
 For both potable and non-potable water, collect two litres of water. Ten standard PHO Laboratory water bottles filled to the 200 mL fill line will satisfy this requirement.
- Complete all sections of the <u>Environmental Microbiology Investigation Requisition</u> form with appropriate information to ensure accurate tracking of sampling and complete all fields as outlined on page 2 of the form.
- Samples should be stored at refrigeration temperature (2 to 8 °C) and shipped to the <u>PHO</u>
 <u>Laboratory</u> as soon as possible as outline in the test directory for <u>Legionella</u> (water).⁵⁴

Considerations

- For potable water systems, e.g., taps or showers, pre-flush samples should be collected to
 represent the water held in the tap or fitting. Do not allow the faucets to run before sampling.
 Turn water on and allow the water to run slowly into the sample bottle to minimize aerosol
 production. Removing the aerator from the faucet can also help to minimize aerosol production.
- Use of Personal Protective Equipment (PPE) should be guided by local board of health internal health and safety guidelines.

Legionella – Laboratory Acceptance Criteria

The accuracy of the test results may be affected by improper collection, handling and/or shipping. *Legionella* samples that do not meet the acceptance criteria may be rejected by the laboratory and a new sample may be re-submitted with a newly completed requisition. The acceptance criteria are outlined in the table below.

Description	Acceptance Reason
Submitter	 Sample must be received from an authorized board of health submitter.
Sample Type	 Swabs and water will be accepted for testing to support investigations of confirmed clinical cases or an outbreak investigation of legionellosis.
Sample Requisition	 The requisition must be completed when it is received at the laboratory including date and location of collection. If a sample is received without a requisition it will not be processed; however, the sample will be held for 24 hours. The requisition must have a unique identifier that matches the identifier on the water collection bottle or swab collection vial.
Sample Collection	 PHO Laboratory water collection bottles and environmental swabs are available from the PHO Laboratory warehouse and can be ordered by calling PHO Laboratory Customer Service Centre at 416-235-6556 or toll free 1-877- 604-4567. Unapproved containers will not be accepted.
Sample Transport	 The sample container must be secure to avoid leaks during transport. The sample temperature must meet the requirements listed below when received at the laboratory: Samples should be shipped with cold packs to maintain a temperature between 2.0 to 8.0 °C when received at the laboratory. Swabs received frozen or greater than room temperature may be subjected to cancellation.

Legionella – Reporting and Interpretation of Environmental Microbiological Test Results

Several factors must be considered when interpreting *Legionella* results, such as sampling points, sampling in relation to water system treatment (use of biocide or thermal shock), conditions to support microbial growth such as temperature and stagnation of water; and any delay or temperature changes in shipping conditions. Collection of water and swabs can be used to describe the colonization of the organism within the system at the time of collection only, recognizing that conditions can change over time. The following table describes the reporting and acceptable limits for *Legionella*.

Testing Type	Reporting Limit	Acceptable Limit
<u>Legionella (water)</u>	Detected / Not detected	Not detected
<u>Legionella (swab)</u>	Detected / Not detected	Not detected

Legionella is ubiquitous to the environment – water, soil and sediment; and can grow in a wide range of temperatures. If Legionella was isolated from an environmental sampling site in connection to an outbreak investigation and is a match to the clinical case, remediation is required and consultation with PHO Laboratory microbiologist is recommended. If Legionella was isolated and is not a match to the clinical case, remediation may be required and PHO Laboratory microbiologist is available for consultation.

Remedial action may be performed immediately if there is strong epidemiological information to suggest potential *Legionella* exposure to a susceptible population. Post remediation samples may be collected and submitted to verify if a previous laboratory confirmed positive source has been effectively treated after remedial action.⁵⁷

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the Public Health Ontario website to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change	
August 8, 2021	 Updated the details of the swab requisition from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition. 	
	 Updated the volume requirement for water testing. 	

Glossary

Specific organism information can be found on the <u>Public Health Agency of Canada website</u>, <u>Pathogen Safety Data Sheet</u>. ⁶⁰ Refer to <u>Appendix D</u>: Water Testing Fact Sheets for interpretation and additional details on Background or Non-target bacteria, Total Coliform, *Escherichia coli*, Heterotrophic Plate count (HPC), *Pseudomonas aeruginosa* and *Staphylococcus aureus*.

Aerobic Colony Count (ACC) – An enumeration of viable aerobic bacteria (bacteria is able grow in an oxygenated or aerobic environment) on microbiological media. This may also be referred to as Standard Plate Count (SPC).

Aerobe – An organism that requires the presence of air or free oxygen to live.

Aerosolization – Dispersion in air of a liquid material or a solution in the form of a fine mist.

Aesthetically unacceptable drinking water – Drinking water that is not clear, colourless or odourless.

Anaerobe – An anaerobic organism lives in the absence of air or free oxygen.

Approval (Directive) – A directive is issued by a Medical Officer of Health that requires a drinking-water system owner/operator to test for one or more health-related parameter. The parameter and maximum concentration are indicated in the directive.

Aseptic technique – A technique where sterile supplies or equipment and clean practices (e.g., hand washing/ sanitation using alcohol based hand rub), are used to prevent bacterial contamination.

Audit sample – A sample taken outside of the routine testing schedule as part of an audit.

Biofilm – A biofilm is a collection of organic and inorganic, living and dead material collected on a surface. In drinking water pipe networks, biofilms can result in reduced disinfectant residuals, increased bacterial levels and other operational issues.

Blowdown – The water that is drained from cooling equipment to remove mineral build-up.

Bottled water – Potable water in commercially sealed containers. In Canada, bottled water is regulated as a food and therefore it must comply with the <u>Food and Drugs Act</u>. ⁶¹ Section 4 of the Act prohibits the sale of foods which contain poisonous or harmful substances and must be for human consumption. A bottled water manufactured from an inadequately treated source would be in contravention of the Food and Drugs Act.

Chain of custody – An unbroken documented chain of events from the time of sample collection until it is processed and completed.

Chromogenic media – Chromogenic media is microbiological media suitable for incubation, differentiation, or selection of different microorganisms as a result of colour production.

Cistern – A reservoir, tank, or container for storing or holding water.

Colony forming units (CFU) – A unit of measurement used to indicate the number of viable microorganisms identified on a specific bacteriological medium for a specific quantity of sample.

Compliance sample – A sample taken to fulfil the requirements of an Ontario drinking water regulation.

Complaint sample (food – no illness) – A food sample submitted for laboratory evaluations that are not associated with an investigation of foodborne illness. Complaints (food with off-odour, off-colour, texture, etc.) usually come from the general public to the local board of health office. The sample could be part of the actual meal, a control sample or a follow-up sample after remedial action.

Control sample (food - illness) – A food sample submitted in conjunction with a suspected foodborne illness incident which were not part of the suspect meal, but which are of the same type and/or prepared and stored under similar conditions as those actually consumed.

Control sample (food – no illness) – A food sample submitted due to a complaint, part of a HACCP investigation or routine surveillance investigation, which are prepared and stored under typical conditions for consumption.

Control sample (water) – A water sample taken prior to a change in a routine or prior to the application of a corrective action.

Cooling tower – A cooling tower is a specialized heat exchanger that allow air and water to meet resulting in the reduction of the water temperature. As a result, evaporation of a small volume of water reduces the temperature of the water being circulated through the tower.

Demineralizer – An instrument, apparatus or chemical used to remove minerals from water.

Designated facility – A facility as described in Ontario Regulation 170/03 Drinking Water Systems i.e., the following facilities: children and youth care, delivery agent care, health care, social care, a children's camp or a school, private school, university, a college of applied arts and technology, or an institution with authority to grant degrees.⁴⁶

Distribution system – The section of a drinking water system that transports drinking water to the end user.

Etiological agent – The organism identified or suspected as causing foodborne or waterborne illness as evidenced through epidemiological information and/or laboratory investigation.

Facultative anaerobe – An organism that can live in the absence as well as in the presence of atmospheric oxygen.

Follow-up samples – Samples taken after remedial action has been taken to determine whether or not the problem has been resolved.

Frozen – Water at or below 0°C or containing obvious ice fragments.

Geometric mean – The geometric mean is a calculation used to average the bacterial levels of E. coli in samples collected from recreational water.

Geometric Mean = ((X1)(X2)(X3).....(Xn))1/n

where X1, X2, etc. represents the individual data points and n is the total number of data points used in the calculation.⁶²

Hazard Analysis Critical Control Point (HACCP) – A systematic preventive approach to food safety from biological, chemical, and physical hazards in production processes that can cause the finished product to be unsafe, and designs measurements to reduce these risks to a safe level.⁶³ For food testing, food samples identified as HACCP on the <u>Food Bacteriology Requisition</u> means, food collected in association with evaluating or monitoring critical control points in food preparation, production and processing.

Heterotrophic Plate Count (HPC) – A microbiological test method used to estimate the number of viable, heterotrophic bacteria in a matrix. It can also be used to measure the changes that occur to bacterial populations during water treatment and distribution. This may also be referred to as Standard Plate Count (SPC).

Hold time – The time from sample collection to time of analysis.

HPPA O. Reg. 318/08 Drinking Water – Drinking water that is subject to the requirements as described in <u>Ontario Regulation 318/08</u> Transitional – Small Drinking Water Systems under the *Health Protection and Promotion Act* until it is repealed.⁶⁴

HPPA O. Reg 319/08 Drinking Water – Drinking water that is subject to the requirements as described in Ontario Regulation 319/08 Small Drinking Water Systems under the *Health Protection and Promotion Act.*⁴⁷

HPPA Regulated Premises – Premises/facilities regulated under HPPA, (i.e., Food premises, recreational water places, etc.) that fall outside the scope of <u>Ontario Regulation 319/08</u> Small Drinking Water Systems.⁴⁷ Refer to the *Examples of Drinking Water that fall outside Ontario Drinking Water regulations* table.

Ice – Potable water which has been cooled to a temperature less than 0 °C.

Illness (as it relates to completion of the Food Bacteriology Requisition) – Foods associated to suspect foodborne illness incidents, are to be submitted using the Food Bacteriology Requisition under the category 'Illness' for reason for test request. These foods can be part of an investigation and/or outbreak; and the samples can be part of the actual meal, control samples or follow-up samples taken after remedial action. All available epidemiological information in conjunction with the foodborne illness incident must also be documented. Epidemiological information includes symptoms, time and date of onset of symptoms, incubation period, and the number of presumed ill and at risk. Documenting the iPHIS case ID number will allow the laboratory to link the food samples to the associated clinical case.

Indicator bacteria (water) – Types of bacteria used to detect and estimate the level of contamination. Some examples include Total Coliforms and *E. coli* (*Escherichia coli*). Indicator bacteria can also be used as process microbial indicators to evaluate the efficacy of a process (e.g., Heterotrophic Plate Count or Total Coliforms for chlorine disinfection), as a fecal indicators (*E. coli*) or index and model organisms which is a group/or species indicative of pathogen presence and behaviour respectively (e.g., *E. coli* as an index for *Salmonella*).⁶⁵

Interactive Voice Response (IVR) – A technology that allows a computer to interact with humans through the use of voice and dual-tone multi-frequency signaling keypad inputs. Private Citizen drinking water test results and the interpretation can be accessed by calling a toll free telephone number and keying in the sample barcode when prompted.

Interfering substances – A substance in a water sample that either prevents the water from passing through the membrane filter during the testing process or post incubation, interferes with the determination of the colony count due to interfering or inhibiting substances on the membrane filter (e.g., a layer of sediment on the membrane filter, etc.).

Legal sample – Any laboratory specimen that has been submitted for testing to the PHO Laboratory with an intact legal seal and for which chain of custody has been maintained.

Legal seal – A tamper evident label attached to a sample or transport container that must be broken before access can be obtained.

Membrane filter – A 47 millimetre diameter cellulose ester membrane with a $0.45 \mu m$ pore size used in the membrane filtration water testing method. A vacuum draws a measured volume of water through the membrane filter enabling the capture of bacteria if present in a sample.

Most Probable Number (MPN) – A method used in microbiology to determine a statistical value that represents the viable microorganisms in a sample through the use of dilution and multiple tube inoculations.

Municipal Drinking-Water System – A municipal drinking water system means a drinking water system or part of a drinking water system, (a) that is owned by a municipality or by a municipal service board established under the <u>Municipal Act, 2001</u> or a city board established under the <u>City of Toronto Act, 2006</u>, (b) that is owned by a corporation established under sections 9, 10 and 11 of the Municipal Act, 2001 in accordance with section 203 of that Act or under sections 7 and 8 of the City of Toronto Act, 2006 in accordance with sections 148 and 154 of that Act, (c) from which a municipality obtains or will obtain water under the terms of a contract between the municipality and the owner of the system, or (d) that is in a prescribed class. 66-67

Non Municipal Drinking Water System – A non-municipal drinking water system means a drinking water system that is not a municipal drinking water system.⁴⁵

Non-regulated drinking water – Drinking water systems that are not captured by specific regulation where the public has access to the water and such systems are subject of an investigation under the

<u>Health Protection and Promotion Act</u>. ⁴⁸ Public health inspectors may collect drinking water samples in concurrence with the following sections of the HPPA: Section 13, Section 41, Subsections (1), (2), (4), and (5). See the table Examples of Drinking Water that fall outside Ontario's Drinking Water regulations.

Official Agency – The following agencies that submit samples for testing to the PHO Laboratory are considered Official Agencies:

- a) Board of health (public health unit): An official health agency established by a group of urban and rural municipalities to provide a more efficient community health program, carried out by full-time, specially qualified staff, or
- b) Ministry of Natural Resources and Forestry Ontario Park: Ontario land dedicated to the people of Ontario and visitors for their inspiration, education, health, recreational enjoyment and other benefits with the intention that these areas shall be managed to maintain their ecological integrity and to leave them unimpaired for future generations.⁶⁸

Order (Provincial Officer's Order) – A mandate made under the <u>Safe Drinking Water Act, 2002</u> or the <u>Health Protection and Promotion Act R.S.O 1990, c.H.7</u> where, in the opinion of the Director or Medical Officer of Health or designate, an imminent drinking water health hazard exists in respect of a regulated drinking water system, the Director/Medical Officer of Health or designate may issue an order to the owner or the operating authority for the system, requiring the owner or operating authority, in accordance with the directions specified in the order. Provincial officer means a person who is designated by the Minister as a provincial officer for the purposes of the Act and the regulations.

Outbreak – The occurrence of disease beyond what is normally expected based on epidemiological and surveillance data.

Outbreak investigation – A sample taken as part of a public health outbreak investigation.

Outbreak number – The number assigned to an outbreak by board of health staff (e.g., public health unit) for tracking purposes.

Overgrown: NDOGHPC - No Data: Overgrown Heterotrophic Plate Count – Crowding and/or confluent microbial growth on standard plate count agar.

Overgrown: NDOGN - No Data Overgrown with non-target – Crowding and/or confluent and/or non-identifiable microbial growth on an agar plate with no evidence of Total Coliform or *E. coli*. The test is considered overgrown with background bacteria, i.e., organisms often found in the environment. This heavy load of background bacteria can interfere with the detection of coliforms and/or *E. coli*. The water may be unsafe to drink.

Overgrown: NDOGT - No Data Overgrown with target – Crowding and /or confluent and/or non-identifiable microbial growth on an agar plate, with evidence of Total Coliforms or Total Coliforms and *E. coli* present. The water is unsafe to drink.

Oxidase – A test performed to differentiate bacteria that produce cytochrome oxidase, an enzyme of the bacterial electron transport chain. When present, the cytochrome c oxidase oxidizes the reagent (tetramethyl-p-phenylenediamine) to (indophenols) and a purple colour end product results. The majority of coliforms are oxidase negative whereas the majority of *Pseudomonas aeruginosa* and *Aeromonas hydrophila* are oxidase positive.

Part of meal – Foods submitted in conjunction with a suspect foodborne illness incident or complaint which is part of the actual meal implicated.

Potable – Fit or suitable for drinking.

Perishable – A type of food with a limited shelf life if it is not refrigerated.

Personal Protective Equipment (PPE) – Protective clothing, helmets, goggles, gloves or other garments or equipment designed to protect the wearer's body from injury or infection.

Prion – A small proteinaceous infectious disease-causing agent that is believed to be the smallest infectious particle; it is neither bacterial nor fungal nor viral and contains no genetic material.

Private Citizen – A person residing in Ontario who:

- a) does not have access to a treated municipal drinking water supply and
- b) draws their drinking water from a protected source within the Province of Ontario for personal use or applies appropriate treatment to water drawn from an unprotected source.

Private Citizen drinking water – Drinking water from a "single household".

Private residence – A dwelling place occupied for an extended period of time by the same persons, if,

- a) the residents have a reasonable expectation of privacy;
- b) food preparation, personal hygiene, and sleeping accommodations are not communal in nature; and
- c) any use of the dwelling place by a resident for a home occupation, trade, business, profession or craft is secondary to the use of the dwelling place as a residence and does not use more than 25 per cent of the indoor floor area.⁶⁹ (e.g., water source can include a private well, cisterns, treated surface water, etc.)

Provincial Officer's order – Refer to Order (Provincial Officer's Order).

Public beach – Any public bathing area owned/operated by a municipality to which the general public has access, and where there is reason to believe that there is recreational use of the water (e.g., beach signage, sectioned off swimming area, water safety/rescue equipment, lifeguard chairs, etc.), which may result in waterborne illness or injury as determined by the local medical officer of health.⁷⁰

Public pool – A structure, basin, chamber or tank containing or intended to contain an artificial body of water for swimming, water sport, water recreation or entertainment, but does not include,

- a) one that is located on a private residential property under the control of the owner or occupant and that is limited to use for swimming or bathing by the owner or occupant, members of their family and their visitors, or
- b) one that is used solely for commercial display and demonstration purposes.⁷¹

Public spa – A hydro-massage pool containing an artificial body of water that is intended primarily for therapeutic or recreational use, that is not drained, cleaned or refilled before use by each individual and that utilizes hydrojet circulation, air induction bubbles, current flow or a combination of them over the majority of the pool area.⁷²

Raw foods – Foods which have NOT been subjected to some form of processing in order to render them RTE (most often cooking) and/or which have been subjected to another process to extend their shelf-life, including but not restricted to the use of heat, chemicals, reduction of pH, reduction of water activity, or special packaging.⁷³

Ready-to-eat foods – The status of the food being ready for immediate consumption at the point of sale. It could be raw or cooked, hot or chilled, and can be consumed without further heat-treatment including re-heating.⁷⁴

Recreational water facility – A recreational water facility can include:

- a) regulated recreational water facilities including public pools and public spas and
- b) non-regulated recreational water facilities including public wading pools, splash pads/spray pads and water slide receiving basins.⁷⁰

Regulated water seal – A tamper evident label attached to the cap of a drinking water sample submitted under a drinking water regulation that must be broken before access can be obtained.

Relinquish – To release a sample ensuring the chain of custody is not compromised.

Resample – A sample taken from the same location, usually in response to an adverse test result or incident.

Routine surveillance (food) – Regular on-going sanitary monitoring programs of food prepared/served in retail food establishments. Samples collected in the absence of foodborne disease.

Sanitizer – A substance or preparation for killing germs, designed for use especially on food-processing equipment.

SDWA O. Reg 170/03 Drinking Water – Drinking water that is subject to the requirements as described in <u>Ontario Regulation 170/03 Drinking Water Systems</u>. ⁴⁶ (e.g., some Municipal, Non-municipal and designated facilities).

Single Household – In the context of a private citizen drinking water, a single household is defined as having a drinking water system that serves one or fewer than six (6) private residences and the water is not available for public consumption.

Spread plate – The technique used to determine the load of heterotrophic bacteria in a water sample. A small volume of sample (0.1 mL) is inoculated on to the surface of an agar plate, spread using aseptic technique and the plate is then incubated according to the test method.

Surface water – Water collecting on the ground or in a stream, river, lake, etc.

Suspected sewage contamination – Water suspected of being contaminated by fecal material.

Teletypewriter (TTY) – A telegraphic apparatus by which signals are sent by striking the letters and symbols of the keyboard of an instrument resembling a typewriter and are received by a similar instrument that automatically prints them in type corresponding to the keys struck.

Test for Potability – Total Coliform and *E. coli* analysis.

Treated water – Water that is subjected to disinfection or other means of removing contaminants.

Turnaround time (TAT) – The number of business days to complete testing based on date of receipt in the laboratory. This turnaround time definition does not apply to tests that are on hold awaiting information to continue testing. Turnaround times posted in this Guide are estimates for routine tests. For some tests or at times of public health emergencies, outbreaks or exigent circumstances, turnaround times may be longer than what is listed.

Unauthorized Submitter: Official Agency – A submitter representing an agency other than a board of health or Ontario Parks under the Ministry of Natural Resources and Forestry.

Unauthorized Submitter: Private Citizen – A person who submits a drinking water sample that does not meet the definition of "Private Citizen" or "Single Household".

Unique identifier – A sample identifier that uniquely distinguishes it from any other sample submitted. For water samples, the sample barcode is used for this purpose.

Untreated – Water that is not subjected to disinfection or other means of removing contaminants.

Water well – A hole drilled, bored or dug into the earth to obtain water.

Acronyms

AOAC Association of Official Analytical Chemists

APHA American Public Health Association

ACC Aerobic Colony Count

ASPHIO Association of Supervisors of Public Health Inspectors of Ontario

a_w Water activity

BOH Board of health

CALA Canadian Association for Laboratory Accreditation

CIPHI Canadian Institute of Public Health inspectors

COC Chain of custody

CFIA Canadian Food Inspection Agency

CFU Colony forming unit

FDA-BAM US Food and Drug Administration Bacteriological Analytical Manual

HACCP Hazard Analysis Critical Control Point

HC Health Canada

HPC Heterotrophic Plate Count

HPPA Health Protection and Promotion Act

iPHIS Integrated Public Health Information System

IU International units

MNRF Ministry of Natural Resources and Forestry

MECP Ministry of the Environment, Conservation and Parks

MOH Medical Officer of Health

MOH-EH Ministry of Health (Ministry of Health, Environmental Health Policy and

Programs, Health Protection and Surveillance Policy and Programs Branch)

NDOGHPC No Data: Overgrown Heterotrophic Plate Count

NDOGN No Data: Overgrown with Non Target

NDOGT No Data: Overgrown with Target

OAHPP Ontario Agency for Health Protection and Promotion

OMAFRA Ontario Ministry of Agriculture, Food, and Rural Affairs

ON-ARC Ontario-Area Recall Coordinator

O. Reg Ontario Regulation

PHI Public health inspector

PHO Public Health Ontario

PHOL Public Health Ontario Laboratory

PHU Public health unit

PPE Personal protective equipment

SCC Standards Council of Canada

SDWA Safe Drinking Water Act

SDWS Small Drinking Water Systems

TAT Turnaround time

UPC Universal Product Code

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change
August 8, 2021	Added hold time definition.
	Updated MOHLTC to MOH-EH.

Appendices

Appendix A: Supplies

Board of health staff should ensure that environmental sample collection and submission supplies are readily available. Collection kits or other PHO Laboratory supplies can be ordered by completing the Requisition for Specimen Containers and Supplies form and faxing it to the PHO Laboratory Toronto fax number listed on the form or to your local PHO Laboratory. Supplies will be shipped within 1 – 5 days of request.

Ordering instructions for supplies that are not listed on the <u>Requisition for Specimen Containers and Supplies</u> form:

ITEM	CATALOGUE #	ORDERING INSTRUCTIONS
Regulated water seal	L-8001 Order directly from local PHO Laborator	
Legal seal	L-8002	Order directly from local PHO Laboratory
Sterile sampling bags with wire closure	300240	Record "Food Bag PT# 300240 B /250" on the comment line of the Requisition for Specimen Containers and Supplies form.
Sewage Collection Kits	TC 125-B	Order directly from local PHO Laboratory
Swab		Record "Swab" on the comment line of the <u>Requisition for Specimen Containers</u> and <u>Supplies</u> form

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change
May 6, 2019	No changes noted.

Appendix B: Laboratory Submission Forms

The most current version of the forms can be accessed from the links below:

Bacteriological Analysis of Drinking Water for Private Citizen, Single Household Only

Bacteriological Analysis of Water - Multiple Sample Requisition for Official Agencies

Environmental Microbiology Investigation Form

Food Bacteriology Requisition Form

Potability of Water Release Form

Regulated water seal (L-8001)

Public Health Ontario Ontario Ontario Ontario Ontario Ontario Ontario Ontario	ublique Ontario	Seal no₌ Numero du Sceau	Note - Place centre of seal over cap so that the bottle cannot be opened without
	D WATER SEAL RUSAGE OFFICIEL		breaking the seal. Ensure that seal is securely fastened.
Initial of person attaching official seal L'initial d'individu apposant le sceau Initial of person receiving specimen	Date		Veuillez vous assurez que le centre du sceau soit placé au dessus du capuchon, et le tout scellé fermement.
Initiales de la personne recevant l'échantillon			10/15 Ontario

Legal seal (L-8002)

Public Health Ontario Ontario Santé publique Ontario Laboratoires de Santé pu		Seal no. Numero du Sceau	Note - Place centre of seal over cap so that the bottle cannot be opened without
	AL SEAL U LÉGAL		breaking the seal. Ensure that seal is securely fastened.
Initial of person attaching official seal L'initial d'individu apposant le sceau	Date		Veuillez vous assurez que le centre du sceau soit placé au dessus du capuchon, et le
Initial of person receiving specimen Initiales de la personne recevant l'échantillon		•	tout scellé fermement. Ontarto TO/15 Ontarto Manufactura Manufact

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the Public Health Ontario website to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change	
August 8, 2021	 Updated the details of the swab requisition from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition. 	

Appendix C: Major Foodborne Diseases: General

Features

Information extracted from World Health Organization (2008). Foodborne Disease Outbreaks: Guidelines for Investigation and Control. Geneva, Switzerland, 2008.⁷⁵

The following tables includes the general features of major foodborne diseases including: etiological agent, incubation period, symptoms, and mode of transmission and associated foods.

BACILLUS CEREUS GASTROENTERITIS

Features	Description
Etiological Agent	Bacillus cereus (Bacterial toxin)
Incubation Period	a) Diarrheal syndrome: 8 – 16 hours b) Emetic syndrome: 1 – 5 hours
Symptoms	a) Diarrheal syndrome: acute diarrhea, nausea and abdominal pain.b) Emetic syndrome: acute nausea, vomiting and abdominal pain and sometimes diarrhea.
Mode of Transmission and Associated Foods	Ingestion of food that has been stored at ambient temperatures after cooking, permitting the growth of bacterial spores and toxin production. Many outbreaks (particularly those of the emetic syndrome) are associated with cooked or fried rice that has been kept at ambient temperature. Foods involved include starchy products such as boiled or fried rice,
	spices, dried foods, milk, dairy products, vegetable dishes, and sauces.

CAMPYLOBACTERIOSIS

Features	Description
Etiological Agent	Campylobacter jejuni and Campylobacter coli
Incubation Period	Typically 2 – 5 days (range 1 – 11 days)

Features	Description
Symptoms	Fever, severe abdominal pain, nausea and diarrhea which can vary from slight to profuse and watery sometimes containing blood or mucus.
Mode of Transmission and Associated Foods	Principally through ingestion of contaminated food. Main food sources are raw milk and raw or undercooked poultry. Spread to other foods by cross-contamination or contamination with untreated water; contact with animals or birds. Other sources of transmission are contact with live animals (pets and farm animals). Person-toperson transmission occurs during the infectious period that ranges from several days to several weeks. Foods involved include raw milk, poultry, beef, pork and drinkingwater

CLOSTRIDIUM BOTULINUM (BOTULISM)

Features	Description
Etiological Agent	Clostridium botulinum
Incubation Period	12 – 36 hours (range several hours to 8 days)
Symptoms	Vomiting, abdominal pain, fatigue, muscle weakness, headache, dizziness, ocular disturbance (blurred or double vision, dilated pupils, unreactive to light), constipation, dry mouth and difficulty in swallowing and speaking, and ultimately paralysis and respiratory or heart failure.
Mode of Transmission and Associated Foods	Ingestion of toxin pre-formed in food. This may occur when raw or under- processed foods are stored in anaerobic conditions that allow growth of the organism. Most outbreaks are due to faulty preservation of food (particularly in homes or cottage industries), e.g. canning, fermentation, curing, smoking, or acid or oil preservation.
	Foods involved include vegetables, condiments (e.g. pepper), fish and fish products, meat and meat products, honey, fruit and vegetable juices. Several outbreaks have occurred as a result of consumption of uneviscerated fish, garlic in oil, and baked potatoes.

CLOSTRIDIUM PERFRINGENS ENTERITIS

Features	Description
Etiological Agent	Clostridium perfringens
Incubation Period	8 – 24 hours
Symptoms	Abdominal pain, diarrhea, rarely vomiting and fever.
Mode of Transmission and Associated Foods	Illness usually caused by cooked meat and poultry dishes subject to time/temperature abuse. Dishes are often left for too long at ambient temperature to cool down before storage, or cooled inadequately. This allows spores that survive the cooking process to germinate and grow, producing large numbers of vegetative cells. If a dish is not reheated sufficiently before consumption, the vegetative cells can cause illness. Foods involved include meat and poultry (boiled, stewed or casseroled).

CRYPTOSPORIDIOSIS

Features	Description
Etiological Agent	Cryptosporidium parvum
Incubation Period	2 – 4 days
Symptoms	Persistent diarrhea, nausea, vomiting and abdominal pain, sometimes accompanied by an influenza-like illness with fever.
Mode of Transmission and Associated Foods	Spread through the fecal-oral route, person-to-person contact or consumption of fecally contaminated food and water, bathing in contaminated pools. Foods involved include raw milk, drinking-water and apple cider.

CYCLOSPORIASIS*

Features	Description
Etiological Agent	Cyclospora cayetanensis
Incubation Period	2 – 11 days
Symptoms	Watery diarrhea (this is the most common symptom), abdominal bloating and gas, fatigue, (tiredness), stomach cramps, loss of appetite, weight loss, mild fever and nausea.
Mode of Transmission and Associated Foods	Cyclosporiasis is unlikely to spread directly between people. This is because the parasite can only infect others once it leaves your body through feces. To be able to spread, the parasite needs to be outside the body for about 7 to 15 days.
	Food can be a source of cyclosporiasis for Canadians when imported from countries where <i>Cyclospora</i> is common. Foods imported to Canada that have been linked to the <i>Cyclospora</i> parasite include: basil, cilantro, raspberries, blackberries, mesclun lettuce, snow and snap peas and pre-packaged salad mix.

^{*}Organism specific reference taken from Government of Canada website⁷⁶

ESCHERICHIA COLI INFECTION

Features	Description
	Escherichia coli
	a) Enteropathogenic <i>E. coli</i> (EPEC)
Etiological Agent	b) Enterotoxigenic <i>E. coli</i> (ETEC) producing a heat-labile (LT) and a heat-stable (ST) enterotoxin
	c) Enteroinvasive E. coli (EIEC)
	d) Enterohaemorrhagic <i>E. coli</i> (EHEC) or verocytotoxin producing <i>E. coli</i> (VTEC), also referred to as Shiga-toxin producing <i>E. coli</i> (STEC), of which the most commonly recognized is <i>E. coli</i> O157.
	a) EPEC: 1 – 6 days; as short as 12 – 36 hours
Incubation Period	b) ETEC: 1 – 3 days; as short as 10–12 hours
	c) EIEC: 1 – 3 days; as short as 10–18 hours

Features	Description
	d) EHEC: 3 – 8 days, median of 4 days
Symptoms	a) EPEC adheres to the mucosa and changes its absorption capacity, causing vomiting, diarrhea, abdominal pain and fever.
	b) ETEC mediates its effects by enterotoxins. Symptoms include diarrhea (ranging from mild to a severe, cholera-like syndrome), abdominal cramps and vomiting, sometimes leading to dehydration and shock.
	c) EIEC causes inflammatory disease of the mucosa and submucosa by invading and multiplying in the epithelial cells of the colon. Symptoms include abdominal pain, vomiting and watery diarrhoea (in $<10\%$ of cases stools may become bloody and contain mucus).
	d) EHEC causes abdominal cramps and watery diarrhea that may also develop into bloody diarrhea (hemorrhagic colitis). Fever and vomiting may also occur.
Mode of Transmission and Associated Foods	(a – c) EPEC, ETEC, EIEC: consumption of food and water contaminated with fecal matter. Time/temperature abuse of such foods increases risk of illness. Up to 25% of infections in infants and young children in developing countries are due to E. coli, in particular ETEC and EPEC ($10 - 20\%$ and $1 - 5\%$ of cases at treatment centres, respectively). ETEC is a major cause of traveller's diarrhea in developing countries.
	d) EHEC is transmitted mainly through consumption of foods such as raw or undercooked ground-meat products and raw milk from infected animals. Fecal contamination of water and other foods, as well as cross-contamination during food preparation, will also lead to infection.
	Foods involved include ground (minced) meat, raw milk, and vegetables. Secondary transmission (person-to-person) may also occur during the period of excretion of the pathogen which is less than a week for adults but up to 3 weeks in one-third of affected children

GIARDIASIS

Features	Description
Etiological Agent	Giardia lamblia
Incubation Period	7 – 10 days (range 4 – 25 days)
Symptoms	Diarrhea (which may be chronic and relapsing), abdominal cramps, fatigue, weight loss, anorexia and nausea. Symptoms may be caused by a protein toxin.
Mode of Transmission and Associated Foods	Infected individuals excrete <i>Giardia</i> cysts in large numbers. Illness is spread by fecal-oral route, person-to-person contact or fecally contaminated food and water. Cysts have been isolated from lettuces and fruits such as strawberries. Infection also associated with drinking-water from surface waters and shallow wells. Foods involved include water, home-canned salmon, fruit and vegetables and noodle salad.

HEPATITIS A

Features	Description
Etiological Agent	Hepatitis A
Incubation Period	25 – 28 days (range 2 – 6 weeks)
Symptoms	Loss of appetite, fever, malaise, abdominal discomfort, nausea and vomiting, followed by symptoms of liver damage (passage of dark urine, pale stools, jaundice).
Mode of Transmission and Associated Foods	Spread by fecal-oral route, primarily person-to-person. Can also be transmitted through food and water as a result of sewage contamination or infected food handlers. Risk of transmission is greatest during the second half of the incubation period until a few days after the appearance of jaundice. Foods involved include shellfish, raw fruit and vegetables, bakery products.

LISTERIOSIS

Features	Description
Etiological Agent	Listeria monocytogenes
Incubation Period	3 – 70 days
Symptoms	Influenza-like symptoms such as fever, headache and occasionally gastrointestinal symptoms.
Mode of Transmission and Associated Foods	A substantial proportion of cases of listeriosis are foodborne. Foods involved include raw milk, soft cheese, meat-based paste, jellied pork tongue, raw vegetables and coleslaw.

SALMONELLOSIS

Features	Description
Etiological Agent	Non-typhoid Salmonella serotypes
Incubation Period	6 – 48 hours, occasionally up to 4 days
Symptoms	The principal symptoms are fever, headache, nausea, vomiting, abdominal pain and diarrhea.
Mode of Transmission and Associated Foods	Main route of transmission is by ingestion of the organisms in food (milk, meat, poultry, eggs) derived from infected food animals. Food can also be contaminated by infected food-handlers, pets and pests, or by cross-contamination as a result of poor hygiene. Contamination of food and water from the feces of an infected animal or person may also occur. Problems caused by initial contamination may be exacerbated by prolonged storage at temperatures at which the organism may grow. Direct person toperson transmission may also occur during the course of the infection. Foods involved include unpasteurized milk, raw eggs, poultry, meat, spices, salads and chocolate.

SALMONELLA TYPHI AND SALMONELLA PARATYPHI TYPES A-C (TYPHOID FEVER, PARATYPHOID FEVER)

Features	Description
Etiological Agent	Salmonella typhi and Salmonella paratyphi types a-c
Incubation Period	10 – 20 days (range 3 days to 8 weeks)
Symptoms	Systemic infections characterized by high fever, abdominal pains, headache, vomiting, diarrhea followed by constipation, rashes and other symptoms of generalized infection.
Mode of Transmission and Associated Foods	Ingestion of food and water contaminated with fecal matter. Foodhandlers may carry the pathogen and be a source of food contamination. Secondary transmission may also occur. Foods involved include prepared foods, dairy products (e.g. raw milk), meat products, shellfish, vegetables, and salads.

SHIGELLOSIS (BACILLARY DYSENTERY)

Features	Description
Etiological Agent	Shigella dysenteriae, S. flexneri, S. boydii, S. sonnei
Incubation Period	1 – 3 days, up to 1 week for S. dysenteriae
Symptoms	Abdominal pain, vomiting, fever, diarrhea ranging from watery (S. sonnei) to dysenteric with bloody stools, mucus and pus (S. dysenteriae and, to a lesser extent S. flexneri and S. boydii).
Mode of Transmission and Associated Foods	Food and water contaminated with fecal matter. Person-to-person transmission through the fecal—oral route is an important mode of transmission. Food can be contaminated by food-handlers with poor personal hygiene or by use of sewage/wastewater for fertilization. Foods involved include uncooked foods that have received extensive handling, such as mixed salads and vegetables, water and raw milk.

STAPHYLOCOCCUS AUREUS INTOXICATION

Features	Description
Etiological Agent	Staphylococcus aureus
Incubation Period	2 – 6 hours
Symptoms	Intoxication, sometimes of abrupt and violent onset. Severe nausea, cramps, vomiting and prostration, sometimes accompanied by diarrhea.
Mode of Transmission and Associated Foods	Consumption of foods containing the toxin. Foods are contaminated by food-handlers. If storage conditions are inadequate, the bacteria may multiply to produce toxin. Intoxication is often associated with cooked food e.g. meat, in which competitive bacteria have been destroyed.
	Foods involved include prepared foods subject to handling in their preparation (ham, chicken and egg salads, cream-filled products, ice cream, cheese).

VIBRIO CHOLERAE O1 AND O139 (CHOLERA)

Features	Description
Etiological Agent	Vibrio cholerae O1 and O139
Incubation Period	1 – 3 days
Symptoms	Profuse watery diarrhea, which can lead to severe dehydration, collapse and death within a few hours unless lost fluid and salt are replaced; abdominal pain and vomiting.
Mode of Transmission and Associated Foods	Food and water contaminated through contact with fecal matter or infected food handlers. Contamination of vegetables may occur through sewage or wastewater used for irrigation. Person-to- person transmission through the fecal-oral route is also an important mode of transmission. Foods involved include seafood, vegetables, cooked rice and ice.

VIBRIO PARAHAEMOLYTICUS GASTROENTERITIS

Features	Description
Etiological Agent	Vibrio parahaemolyticus
Incubation Period	9–25 hours, up to 3 days
Symptoms	Profuse watery diarrhea, abdominal pain, vomiting, and fever. A dysenteric syndrome has been reported from some countries, particularly Japan.
Mode of Transmission and Associated Foods	Associated with consumption of raw or undercooked fish and fishery products or cooked foods subject to cross contamination from raw fish.

VIBRIO VULNIFICUS INFECTION

Features	Description
Etiological Agent	Vibrio vulnificus
Incubation Period	12 hours – 3 days
Symptoms	Profuse diarrhea with blood in stools. Organism is associated with wound infections and septicaemia may originate from the gastrointestinal tract or traumatized epithelial surfaces.
Mode of Transmission and Associated Foods	All known cases are associated with seafood, particularly raw oysters.

YERSINIOSIS

Features	Description
Etiological Agent	Yersinia enterocolitica
Incubation Period	24 – 36 hours (range 1 – 11 days)
Symptoms	Abdominal pain, diarrhea, mild fever, sometimes vomiting.

Features	Description
Mode of Transmission and Associated Foods	Illness is transmitted through consumption of pork products (tongue, tonsils, gut), cured or uncured, as well as milk and milk products.

VIRAL GASTROENTERITIS

Features	Description
	Viral gastroenteritis
Etiological Agent	Many different viruses can cause viral gastroenteritis, including adenoviruses, coronaviruses, rotaviruses, parvoviruses, caliciviruses and astroviruses. Those viruses most commonly associated with foodborne outbreaks are norovirus and hepatitis A.
Incubation Period	15 – 50 hours
Symptoms	Diarrhea and vomiting, which is often severe and projectile with sudden onset.
Mode of Transmission and Associated Foods	Gastroenteritis viruses usually spread by fecal—oral route. Food and drinking-water may be contaminated either at source when exposed to sewage/wastewater in the environment or used for irrigation, or by an infected food-handler. Filter-feeding shellfish most common food contaminated at source, but a wide range of different cooked and uncooked foods have been implicated in secondary contamination by food-handlers.

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the Public Health Ontario website to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change
May 6, 2019	No changes noted

Appendix D: Water Testing Fact Sheets

The following tables includes the general information for water testing including information about the organism, distinguishing features, where the organism came from, potential health impacts and impact/importance. Reporting and acceptable limits are also presented for each organism.

BACKGROUND BACTERIA / NON TARGET BACTERIA / OVERGROWN

Features	Description
What is it?	Background bacteria are bacteria known as non-target organisms in drinking water testing.
Distinguishing features	These bacteria do not ferment lactose when grown on culture media designed for coliform detection unlike the typical the indicator organisms, coliforms and <i>E. coli</i> .
Where do they come from?	Background bacteria are found in the environment (e.g., soil) and may include potentially "opportunistic" bacterial pathogens most of which are non- coliforms (e.g., <i>Pseudomonas</i> , <i>Aeromonas</i> , etc.) and are not often found in drinking water which has been properly disinfected. ⁷⁸
Potential health impacts	There is no evidence of an association of any of these organisms with gastrointestinal infection through ingestion of drinking-water in the general population. ⁷⁷
	When water samples are heavily contaminated (i.e., overgrown where the plate is crowded and/or confluent and/or contains non-identifiable microbial growth) with background bacteria these organisms can interferes with the detection of coliforms or <i>Escherichia coli</i> that may be present in the sample and in some cases can release inhibitory substances (bacteriocins). ⁷⁸⁻⁷⁹
	Drinking Water
What does it mean?	There are two overgrown conditions:
	 NDOGN - No Data Overgrown with non-target: The crowding and /or confluent and/or non-identifiable microbial growth (i.e. non-target organisms) on a plate, and there is no evidence of coliforms and / or <i>E. coli</i>. The water may be unsafe to drink. NDOGT - No Data Overgrown with target: The crowding and /or confluent and/or non-identifiable microbial growth (i.e. non-target organisms) on a plate, and

Features	Description
	there is evidence of coliforms and / or <i>E. coli</i> present. The water is unsafe to drink.
	If <i>Aeromonas spp</i> . or <i>Pseudomonas aeruginosa</i> is detected in a regulated drinking water sample, it is considered an adverse result for the purpose of section 18 of the <u>Safe Drinking Water Act</u> . ⁴⁵
	Recreational Water Facilities
	An overgrown result is indicative of poor water quality.

REPORTING AND ACCEPTABLE LIMITS: OFFICIAL AGENCY POTABLE WATER / ICE – TREATED / RECREATIONAL WATER FACILITIES

Testing Type	Reporting Limit	Acceptable Limit
Background or Non-target bacteria (CFU / 100 mL)	d or Non-target 0 to *NDOGN /NDOGT mL)	At the discretion of the MOH. Only reported on Regulated drinking water samples submitted under a "Provincial Officer's Order" or "Approval (Directive)". An overgrown result is considered an adverse condition under Ontario drinking
		water regulations.

^{*}NDOGN - No Data Overgrown with non-target / NDOGT- No Data Overgrown with target

CAMPYLOBACTER SPECIES

Features	
Organism	Campylobacter species
What is it?	Campylobacter is a genus of Gram-negative bacteria.
Distinguishing features	Campylobacter are Gram-negative, oxidase and catalase positive and motile. They are microaerophilic (prefer environments with low oxygen) and are best cultured at 42 °C. They generally appear comma or s-shaped, but when exposed to atmospheric oxygen, the pathogenic strain, <i>C. jejuni</i> , is able to change into a coccal form. ⁸⁰
Where do they come from?	Campylobacter jejuni can be isolated from animals with poultry, wild birds and cattle being significant reservoirs. ⁷⁷

Features	
Potential health impacts	Refer to Appendix C.
	Water
	Analysis of water samples for <i>Campylobacter</i> should only be undertaken when there has been clinical confirmation of an outbreak and a water source has been incriminated epidemiologically.
	Potable water
What does it mean?	Drinking-water supplies contaminated with <i>Campylobacter</i> spp. are usually associated with unchlorinated or inadequately chlorinated surface water supplies. Since <i>Campylobacter</i> spp. are fecally borne pathogens and are not particularly resistant to disinfection, <i>E. coli</i> is an appropriate indicator for the presence/absence of <i>Campylobacter</i> spp. in drinking-water supplies. ⁷⁷
	Recreational water
	The occurrence of <i>Campylobacter</i> in surface waters has proved to be strongly dependent on rainfall, water temperature and the presence of waterfowl. ⁸¹

REPORTING AND ACCEPTABLE LIMITS: CAMPLYOBACTER SPECIES

Sample Type	Reporting Limit (per volume of water tested)	Acceptable Limit (per volume of water tested)
Water	Not Detected / Detected	Not detected

Coliform

Features	Description
What is it?	Coliforms are non-specific indicators of bacterial contamination.
Distinguishing features	They are a group of Gram-negative, non spore-forming, rod shaped bacteria capable of fermenting lactose with the production of acid and gas, within 48 hours, when grown in a lactose-containing culture medium incubated at 36+/-1.0 °C. Some coliforms are slow lactose fermenters that can produce the enzyme β -D galactopyranosidase that hydrolyzes chromogenic substrates, such as ortho-nitrophenyl- β -D-

Features	Description
	galactopyranoside (ONPG) or chlorophenol red- β -D-galactopyranoside (CPRG). This reaction generally occurs within 24 +/- 2 h or less at 35 °C. <i>Escherichia coli</i> and thermotolerant coliforms are a subset of the Total Coliform group.
	There are some lactose-fermenting strains of species within the genus <i>Aeromonas</i> which produce a coliform-like reaction on media used for the detection of bacterial indicators of contamination. <i>Aeromonas</i> species are oxidase-positive and generally fail to produce gas during lactose fermentation in liquid media. They do not, therefore, conform to the classical definition of a coliform. Regardless, they are often included in the Total Coliform count because the growth medium fails to differentiate this organism from the true coliforms.
Where do they come from?	Coliforms occur naturally in soil and decaying vegetation, but may also be associated with human or animal fecal contamination, so the presence of coliforms does not necessarily provide specific evidence of fecal contamination. ⁸²
Potential health impacts	The presence of coliforms indicates the water source is not secure increasing an individual's health risk if the water is consumed.
	Potable water
What does it mean?	Coliforms may be indicative of a contaminated water supply and/or inadequate treatment. The presence of Total Coliforms in drinking water distribution systems and stored water supplies may reveal regrowth and possible biofilm formation or contamination from foreign material, including soil or plants.
	In private drinking water systems, follow-up actions should be taken to investigate the cause of the positive results. ⁸³ Given the vulnerability of well water to external influences, such as human, environmental and agricultural impacts, it is important to test water frequently.
	Bottled water
	Although rarely found in bottled water, coliforms are used as an indicator of insanitation or possible contamination. ⁸²
	The Food and Drug Regulations specify that prepackaged water represented as mineral water or spring water and prepackaged ice, shall not contain any coliform bacteria ⁸⁴

Features	Description
	Recreational water facilities
	Total Coliform testing on its own may be insufficient to judge the microbiological quality of whirlpool water; it is advisable to consider the outcome of all the test results for a sampling session. ⁷⁸

REPORTING AND ACCEPTABLE LIMITS: COLIFORMS

Sample Type	Reporting Limit (CFU/100 mL)	Acceptable Limit (CFU / 100 mL)
Official Agency Potable Water including bottled water / Ice – Treated	0 to > 80, *NDOGT	0
Private Citizen – Private Drinking Water Systems	0 to > 80, *NDOGT	≤ 5
Official Agency - Recreational Water Facilities, Public Pools/Spas	< 2 to > 160, *NDOGT	< 2

^{*}NDOGT- No Data Overgrown with target

ESCHERICHIA COLI (E. COLI)

Features	Description	
What is it?	E. coli is an indicator of fecal contamination in water.	
Distinguishing features	<i>E. coli</i> are Gram-negative, facultative anaerobic, non spore-forming, lactose- fermenting, rod-shaped bacteria which are oxidase negative. When using chromogenic media, <i>E. coli</i> are defined as coliform bacteria with the enzyme β-glucuronidase. This enzyme splits the substrate 5-bromo-4-chloro-3- indolyl-β-D-glucuronide (i.e., BCIG) within 24-72 h at 36 °C and hydrolyzes 4- methylumbelliferyl-β-D-glucuronide (MUG) to yield a fluorogenic product that is detectable under long-wave (366 nm) UV light within 4 to 24 hours. Approximately 97% of <i>E. coli</i> strains can produce the enzyme β-D- glucuronidase. 85	

Features	Description
Where do they come from?	The majority of <i>E. coli</i> strains are not pathogenic and live in the intestinal tract of animals, including humans, at high densities $(10^7 \text{ to } 10^8 \text{ cells per gram of stool})$.
Potential health impacts	Some <i>E. coli</i> strains can cause diarrhea when contaminated water are consumed, or through contact with contaminated animals or human feces. Refer to Appendix C - Major Foodborne Diseases: General Features.
	The presence of <i>Escherichia coli</i> in a water source may indicate fecal contamination from a human or animal source and thus the possible presence of fecal pathogenic bacteria, viruses and protozoa which can cause gastrointestinal illness and in some instances potential health complications (e.g., hemolytic uremic syndrome from <i>E. coli</i> O157:H7).
	Although <i>E. coli</i> is a good indicator of the presence of enteric bacterial pathogens, it is important to realize enteric viruses and protozoa are more resistant to disinfection than bacteria, the absence of <i>E. coli</i> does not necessarily mean these organisms are not present. Even with this limitation, <i>E. coli</i> continues to be one of the targets of choice to monitor water quality by the World Health Organization and regulatory bodies. ^{75,83}
	Potable / bottled water and treated ice
What does it mean?	E. coli should not be present in drinking water, including bottled/prepackaged water and ice. 51,82,84,87
	Recreational water
	Studies have shown a strong correlation between the concentration of <i>E. coli</i> in fresh waters and the risk of gastrointestinal illness among swimmers, but enteric pathogens, especially enteric pathogenic viruses and protozoa, may be present in water even when <i>E. coli</i> is absent. ^{83,88,89} A guideline value of ≤ 200 <i>E. coli</i> CFU/100 mL is used in Ontario and is based on a geometric mean concentration of a minimum of five public beach samples. ⁹⁰ Refer to the Ministry of Health and Ministry of Long-Term Care Operational Approaches for Recreational Water Guideline document for direction when test results exceed the recreational water level for <i>E. coli</i> .
	Sewage

Features	Description	
	<i>E. coli</i> levels in water >10 ⁵ CFU per 100 mL are indicative of fecal contamination. ⁹¹	

REPORTING AND ACCEPTABLE LIMITS: ESCHERICHIA COLI (E. COLI)

Sample Type	Reporting Limit (CFU/100 mL)	Acceptable Limit (CFU / 100 mL)
Official Agency Potable Water including bottled water / Ice – Treated	0 to > 80, *NDOGT	0
Private Citizen – Private Drinking Water Systems	0 to > 80, *NDOGT	0
Official Agency - Public Beach Water	< 10 to > 1.0 x 10 ³	≤ 200
Official Agency – Suspected Sewage Contamination	< 1.0 x 10 ³ to > 1.0 x 10 ⁶	Not Applicable
Official Agency - Recreational Water Facilities, Public Pools/Spas	< 2 to > 160, *NDOGT	< 2

^{*}NDOGT- No Data Overgrown with target

HETEROTROPHIC PLATE COUNT (HPC)

Features	Description
What is it?	The Heterotrophic Plate Count (HPC) is a process of estimating the number of viable, heterotrophic bacteria in a sample, i.e., bacteria that use organic (carbon-containing) compounds as a source of energy and carbon, and measure the changes that occur to bacterial populations during water treatment and distribution.
Distinguishing features	The non-selective medium used for the HPC test includes carbohydrate, carbon and nitrogen sources that can sustain the growth of a wide variety of organisms in the presence of oxygen. ⁹²

Features	Description	
Where do they come from?	These organisms are ubiquitous and can be found in the environment.	
Potential health impacts	For healthy individuals, drinking water that contains heterotrophic bacteria does not necessarily pose a public health risk. ⁹³	
	Potable water / ice Water with an HPC result ≤ 500/mL is deemed to have detectable	
	residual disinfectant so an HPC > 500 per 1 mL in a distribution system or other treated water is indicative of poor water quality and possible problems with the disinfection procedure. 94 The test has little value as an indicator of pathogen presence but can be useful in operational monitoring as a treatment and disinfectant indicator. In addition, HPC measurement can be used in assessing the cleanliness and integrity of distribution systems and the presence of biofilms. 94	
What does it mean?	Bottled water	
	Water in sealed containers, other than water represented as mineral water or spring water, shall not contain more than 100 total aerobic bacteria (i.e., heterotrophic bacteria) per millilitre as determined by an HPC test. 95 When an HPC test shows unacceptable results, an investigation should be carried out to determine the cause, and that corrective action taken. 96	
	Recreational water facility	
	For public swimming pools, not more than 15 % of samples collected during any 30 day period should have an HPC count of 200 CFU/mL. ⁷⁸	

REPORTING AND ACCEPTABLE LIMITS: HETEROTROPHIC PLATE COUNT (HPC)

Sample Type	Reporting Limit (CFU per 1 mL)	Acceptable Limit (CFU per 1 mL)
Official Agency Potable Water including bottled water / Ice - Treated	< 10 to > 3.0 x10 ³ , * NDOGHPC	< 500
Official Agency - Recreational Water Facilities, Public Pools/Spas	< 10 to > 3.0 x10 ³ , * NDOGHPC	Refer to Appendix D: Water Testing Fact Sheets

^{*}NDOGHPC: No Data: Overgrown Heterotrophic Plate Count

LEGIONELLA SPECIES

Features	Description
What is it?	Legionella is a bacterium. L. pneumophila is the most common cause of legionellosis which can present as two forms of respiratory illness in humans: Legionnaires' disease causing pneumonia and the less severe Pontiac fever that causes a self-limited flu-like illness. ⁹⁷
Distinguishing features	Legionella are gram-negative, aerobic, rod-shaped bacteria that are fastidious and require specific nutrients in order to grow on bacteriologic medium. 97
Where do they come from?	Water is the major reservoir for legionellae. Stagnant water, warm water (between 20°C and 50°C) and biofilm, scale and sediment are all favourable environments for this organism. ⁹⁸
Potential health impacts	They are commonly found in natural water sources (e.g., freshwater lakes, rivers, creeks and hot springs) and can be found in poorly maintained artificial water systems (i.e., cooling towers, air conditioning evaporative condensers, hot and cold water systems and whirlpool spas). Legionellae survive in these aquatic environments as intracellular parasites of free-living protozoa and within biofilms. ⁹⁹
What does it mean?	If <i>Legionella</i> was isolated from an environmental sampling site in connection to an outbreak investigation and matches the serotype of the clinical case, remediation is required and consultation with PHO Laboratory microbiologist is recommended. If <i>Legionella</i> was isolated and is not a match to the clinical case, remediation may be required and PHO Laboratory microbiologist is available for consultation. Refer to PHO Legionellosis website for supporting documentation. 100

REPORTING AND ACCEPTABLE LIMITS: LEGIONELLA SPECIES

Sample Type	Reporting Limit	Acceptable Limit
Water	Detected / Not detected	Not detected

PSEUDOMONAS AERUGINOSA

Features	Description
What is it?	Pseudomonas aeruginosa is an opportunistic bacterial pathogen that can cause skin infections when water is contaminated with significant levels of the bacteria.
Distinguishing features	Pseudomonas aeruginosa is a gram negative, non-fermentative, oxidase positive, rod shaped bacteria, which grows well at 42 °C, is actively motile, and under appropriate incubation conditions produces a fluorescent pigment and often a visible pigment.
	Pseudomonas aeruginosa can be found in feces, soil, natural waters such as lakes and rivers in concentrations of 10/100 mL to > 1,000/100 mL and sewage. They can multiply in water environments and on the surface of suitable organic materials in contact with water (e.g., in drinking water distribution systems and/or biofilms). 101,77,102
Where do they come from?	Pseudomonas aeruginosa can be found in swimming pools, waterslides and is frequently isolated from whirlpools and hot tubs when the organism is shed by infected human or when the organisms is transferred from a bather's feet and hands to the water. It thrives in warm temperatures, and aeration of the water can cause sloughing off of skin cells which also enhances its growth; these cells can act as a nutrient source and can reduce the overall residual disinfectant level as the organic load increases. 103-104
	The risk of colonization from ingesting <i>Pseudomonas aeruginosa</i> in drinking water is low. ¹⁰¹
Potential health impacts	There are two routes of transmission of <i>Pseudomonas aeruginosa</i> in recreational water facilities:
	 Skin exposure, a primary cause of ear, eye and skin infections (e.g., folliculitis) among swimmers⁷⁷ and Lung exposure from inhaling aerosols¹⁰¹
	One study reported the infective dose for healthy individuals is greater than 1,000 organisms per $\rm ml.^{105}$
	Potable / bottled water
What does it mean?	High numbers of <i>Pseudomonas aeruginosa</i> in potable water, particularly packaged water, can result in off- taste, odour and turbidity. Ensuring adequate disinfection and taking steps to control

Features	Description	
	biofilm growth can reduce the growth of <i>Pseudomonas aeruginosa</i> . At source and during marketing, a natural mineral water, spring water or other bottled waters offered for sale is are to be free from <i>Pseudomonas aeruginosa</i> in any 250 ml sample. ^{77,106} If <i>Pseudomonas aeruginosa</i> is detected in a regulated drinking water sample, it is considered an adverse result for the purpose of section 18 of the <u>Safe Drinking Water Act</u> since it is a strong indicator the drinking water system lacks appropriate disinfection or a system malfunction has occurred (i.e., backflow or pipe breakage). ^{45,107}	
	Recreational water facilities	
	<i>Pseudomonas aeruginosa</i> is the primary indicator of disinfection efficacy, while Total Coliforms, heterotrophic plate count, and staphylococci are considered supporting indicators of water quality for these water sources, and is frequently isolated from whirlpool water that is coliform-negative. ⁷⁸	
	Well operated pools should not normally contain <i>Pseudomonas</i>	

Well operated pools should not normally contain *Pseudomonas aeruginosa*. If the count is over 10 CFU of *Pseudomonas aeruginosa* per 100 mL, repeat testing should be undertaken. Where repeated samples contain unacceptable amounts of *Pseudomonas aeruginosa* filtration and disinfection processes should be examined to determine whether there are areas within the pool circulation where the organism is able to multiply.¹⁰⁸

REPORTING AND ACCEPTABLE LIMITS: PSEUDOMONAS AERUGINOSA

Sample Type	Reporting Limit (CFU per100 mL)	Acceptable Limit (CFU per 100 mL)
Official Agency Potable Water including bottled water	0 to > 100 < 2 to > 200 *	0
Official Agency - Recreational Water Facilities, Public Pools/Spas	0 to > 100 < 2 to > 200 *	< 10

^{*} When 50 mL is analyzed

SALMONELLA SPECIES

Features	Description
What is it?	Salmonella is a bacterium.
Distinguishing features	Salmonella are rod-shaped gram-negative bacteria of the family Enterobacteriaceae.
Where do they come from?	Salmonella are carried in the intestinal tract of domestic and wild animals such as cattle, swine, poultry, wild birds, and pets (particularly reptiles) as well as flies. For Salmonella typhi, humans are the only known host. 109
	Since <i>Salmonella</i> cells can be injured or present in low numbers in aqueous environments, large sample volumes are necessary to allow for a concentration procedure as part of the analysis. ¹¹⁰
Potential health impacts	Refer to Appendix C.
	Water
What does it mean?	Salmonella sp. is a pathogen and should not be present in potable water. It is relatively sensitive to disinfection so Escherichia coli is a generally reliable index for Salmonella spp. in drinking-water supplies. ⁷⁷

REPORTING AND ACCEPTABLE LIMITS: SALMONELLA SPECIES

Sample Type	Reporting Limit	Acceptable Limit
Water	Not Detected / Detected	Not detected

STAPHYLOCOCCUS AUREUS

Features	Description
What is it?	Staphylococcus aureus is an opportunistic bacterial pathogen.
Distinguishing features	Gram-positive, non-motile, non-spore-forming coccus, catalase negative facultative anaerobe of the family Micrococcaceae that produces a wide array of toxins and can coagulate rabbit plasma.

Features	Description
Where do they come from?	Staphylococcus aureus are typically coagulase-positive which differentiates them from non-hazardous coagulase-negative strains that are found on skin.
Where do they come from?	Healthy individuals can <i>carry Staphylococcus aureus</i> in their nose, on their skin and in their feces and shed them when they immerse themselves in water. ⁷⁷
	Refer to Appendix C.
Potential health impacts	<i>S. aureus</i> is a major pathogen responsible for purulent infections of the skin (e.g., folliculitis), eyes and ears and may not become apparent until 48 hours after contact with the contaminated water. ⁷⁷
	Potable / bottled water
	There is no evidence of disease transmission of <i>S. aureus</i> through secondary to the consumption of drinking water. Since fecal material is not usually the source of <i>S. aureus</i> , <i>E. coli</i> is not a suitable indicator for the presence of <i>S. aureus</i> in drinking-water. ⁷⁷
What does it mean?	If <i>S. aureus</i> is detected in a regulated drinking water sample, it is considered an adverse result for the purpose of section 18 of the <u>Safe Drinking Water Act</u> since it is a strong indicator the drinking water system lacks appropriate disinfection or a system malfunction has occurred (e.g., backflow or pipe breakage). ^{45,107}
	Recreational water facilities
	The presence of <i>S. aureus</i> in swimming pools and spas is a good index of the health hazard associated with a high bather density, and accounts for a large percentage of swimming-pool-associated illnesses.
	Whenever swimming pool samples are examined for total staphylococci or <i>Staphylococcus aureus</i> , not more than 50 organisms/100 mL should be present. Maintaining free chlorine levels greater than 1 mg/L or equivalent disinfection efficiency should provide adequate inactivation of <i>S. aureus</i> in swimming pools. ^{78,77,111}

REPORTING AND ACCEPTABLE LIMITS: STAPHYLOCOCCUS AUREUS

Sample Type	Reporting Limit (CFU per100 mL)	Acceptable Limit (CFU per 100 mL)
Official Agency Potable Water including bottled water	0 to >100 <2 to >200 *	0
Official Agency - Recreational Water Facilities, Public Pools/Spas: Staphylococci (Presumptive or <i>S. aureus</i>)	0 to >100 <2 to >200 *	<50 **

^{*} When 50 mL is analyzed

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of Implementation	Description and Change
May 6, 2019	 Changed acceptable level of E. coli for public beaches from 100 CFU/100 mL to 200 CFU/100 mL, as per Operational Approaches for
	Recreational Water Guideline, 2018.

^{**} For total staphylococci or Staphylococcus aureus

Appendix E: Ontario Public Health Standards and Supporting Legislation

The Ontario Public Health Standards are currently being updated as the Standards for Public Health Programs and Services (SPHPS). Associated protocols, guidelines and indicators are also undergoing review. Further to release of the updated SPHPS by the Ministry of Health and Ministry of Long-Term Care, these references will be updated and the implications of any changes will be considered throughout the Guide.

ONTARIO PUBLIC HEALTH STANDARDS

 Ontario. Ministry of Health; Ministry of Long-Term Care. Ontario public health standards: requirements for programs, services, and accountability (standards) [Internet]. Toronto, ON: Queen's Printer for Ontario; 2018 [cited 2021 Mar 03]. Available from: health.gov.on.ca/en/pro/programs/publichealth/oph_standards/default.aspx

ONTARIO PUBLIC HEALTH PROTOCOLS

- Ontario. Ministry of Health and Long-Term Care. Food safety protocol, 2019. Toronto, ON:
 Queen's Printer for Ontario; 2019. Available from:
 health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/protocols_guidelines/Food_d_Safety_Protocol_2019_en.pdf
- Ontario. Ministry of Health and Long-Term Care. Safe drinking water and fluoride monitoring protocol, 2019. Toronto, ON: Queen's Printer for Ontario; 2019. Available from:
 health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/protocols_guidelines/Safe
 e Water Fluoride Protocol 2019 en.pdf
- Ontario. Ministry of Health and Long-Term Care. Recreational water protocol, 2019. Toronto,
 ON: Queen's Printer for Ontario; 2019. Available from:
 health.gov.on.ca/en/pro/programs/publichealth/oph standards/docs/protocols guidelines/Rec
 reational Water%20Protocol 2019 en.pdf

ONTARIO PUBLIC HEALTH STANDARDS GUIDELINE DOCUMENTS

- Ontario. Ministry of Health and Long-Term Care. Operational approaches for food safety guideline, 2018. Toronto, ON: Queen's Printer for Ontario; 2019. Available from: health.gov.on.ca/en/pro/programs/publichealth/oph standards/docs/protocols guidelines/Operational Approaches For Food Safety Guideline 2019 en.pdf
- Ontario. Ministry of Health and Long-Term Care. Operational approaches for recreational water guideline, 2018. Toronto, ON: Queen's Printer for Ontario; 2018. Available from: health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/protocols_guidelines/Operational_Approaches_to_Rec_Water_Guideline_2018_en.pdf

- Ontario. Ministry of Health and Long-Term Care. Food safety: a guide for Ontario's food handlers. Toronto, ON: Queen's Printer for Ontario; 2018. Available from: health.gov.on.ca/en/pro/programs/publichealth/enviro/docs/training_manual.pdf
- Ontario. Ministry of Health and Long-Term Care. Response to adverse drinking water quality incidents [Internet]. Toronto, ON: Queen's Printer for Ontario; 2009 [cited 2021 Mar 03].
 Available from: https://collections.ola.org/mon/23007/293938.pdf
- Ontario. Ministry of Health and Long-Term Care. Small drinking water systems risk assessment guideline, 2018. Toronto, ON: Queen's Printer for Ontario; 2018. Available from:
 <u>health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/protocols_guidelines/Sm_all_Drinking_Water_Systems_Risk_Assessment_Guideline_2018_en.pdf</u>

HEALTH PROTECTION AND PROMOTION ACT AND REGULATIONS

Please note: These links are provided for your convenience. To access official copies of Ontario's statutes and regulations, please visit: ontario.ca/laws

- Health Protection and Promotion Act, RSO 1990, c H.7. Available from: ontario.ca/laws/statute/90h07
 - Small Drinking Water Systems, O Reg 319/08. Available from: e-laws.gov.on.ca/html/regs/english/elaws regs 080319 e.htm
 - Public Spas, O Reg 428/05. Available from: e-laws.gov.on.ca/html/regs/english/ elaws regs 050428 e.htm
 - Food Premises, O Reg. 493/17. Available from: ontario.ca/laws/regulation/170493
 - Public Pools, RRO 1990, Reg 565. Available from: e-laws.gov.on.ca/html/regs/english/elaws regs 900565 e.htm

SAFE DRINKING WATER ACT AND REGULATIONS

Please note: These links are provided for your convenience. To access official copies of Ontario's statutes and regulations, please visit: ontario.ca/laws

- Safe Drinking Water Act, SO 2002, c 32. Available from: <u>e-laws.gov.on.ca/html/statutes/english/elaws_statutes_02s32_e.htm</u>
- Ontario Drinking Water Quality Standards, O Reg 169/03. Available from: e-laws.gov.on.ca/ html/regs/english/elaws regs 030169 e.htm
- Drinking Water Systems, O Reg 170/03. Available from: e-laws.gov.on.ca/html/regs/english/ elaws_regs_030170_e.htm

- Definitions of Words and Expressions Used in the Act, O Reg 171/03. Available from: e-laws.gov.on.ca/html/regs/english/elaws_regs_030171_e.htm
- Schools, Private Schools and Day Nurseries, O Reg 243/07. Available from: <u>e-laws.gov.on.ca/html/regs/english/elaws_regs_070243_e.htm</u>
- Drinking Water Testing Services, O Reg 248/03. Available from: e-laws.gov.on.ca/html/regs/english/elaws regs 030248 e.htm

OTHER RELEVANT LEGISLATION

Please note: These links are provided for your convenience. To access official copies of Ontario's statutes and regulations: ontario.ca/laws

- Environmental Protection Act, RSO 1990, c E.19. Available from: e-laws.gov.on.ca/html/ statutes/english/elaws statutes 90e19 e.htm#BK7
- Food Safety and Quality Act, 2001, S.O. 2001, c. 20. Available from: ontario.ca/laws/statute/01f20
- Ontario Water Resources Act, 1990. RSO 1990, c 1990. c O.40. Available from: ontario.ca/laws/statute/90o40
 - Wells, RRO 1990, Reg 903. Available from: e-laws.gov.on.ca/html/regs/english/ elaws regs 900903 e.htm
- Provincial Parks and Conservation Reserves Act, SO 2006, c 12. Available from: elaws.gov.on.ca/html/ statutes/english/elaws_statutes_06p12_e.htm

FOODBORNE ILLNESS OUTBREAK GUIDANCE DOCUMENTS

- Health Canada; Public Health Agency of Canada; Canadian Food Inspection Agency. Weight of
 evidence: factors to consider for appropriate and timely action in a foodborne illness outbreak
 investigation. Ottawa, ON: Her Majesty the Queen in Right of Canada, represented by the
 Minister of Health; 2011. Available from: https://www.canada.ca/content/dam/hc-sc/migration/hc-sc/fn-an/alt_formats/pdf/pubs/securit/2011-food-illness-outbreak-eclosion-malad-ailments-eng.pdf
- Ontario. Ministry of Health. Ontario's Foodborne Illness Outbreak Response Protocol (ON-FIORP) 2020 [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [cited 2021 Mar 03].
 Available from: health.gov.on.ca/en/pro/programs/publichealth/enviro/docs/fiorp_protocol.pdf
- Public Health Agency of Canada. Canada's Foodborne Illness Outbreak Response Protocol (FIORP): a guide to multi-jurisdictional enteric outbreak response. Ottawa, ON: Her Majesty the Queen in Right of Canada, as represented by the Minister of Health; 2017. Available from:

www.canada.ca/content/dam/phac-aspc/documents/services/publications/health-risks-safety/64-02-17-1879-FIORP-2015-EN-04.pdf

SUPPORTING DOCUMENTATION

- American Public Health Association. Compendium of methods for the microbiological examination of foods. 5th ed. Washington, DC: American Public Health Association; 2015.
- Centers for Disease Control and Prevention. Sampling procedure and potential sampling sites
 [Internet]. Atlanta, GA: Centers for Disease Control and Prevention; 2015 [cited 2021 Mar 04].
 Available from: cdc.gov/legionella/downloads/cdc-sampling-procedure.pdf
- Health Canada, Health Products and Food Branch (HPFB). The compendium of analytical methods [Internet]. Evergreen ed. Ottawa, ON: Government of Canada; 1989 [modified 2018 May 09; cited 2021 Mar 04]. Available from: https://www.canada.ca/en/health-canada/services/food-nutrition/research-programs-analytical-methods/analytical-methods/compendium-methods.html
- Heymann DL, ed. Control of communicable diseases manual. 20th ed. Washington, DC: American Public Health Association; 2015.
- Ontario Agency for Health Protection and Promotion (Public Health Ontario). Home canning: literature review. Toronto, ON: Queen's Printer for Ontario; 2014. Available from: https://www.publichealthontario.ca/-/media/documents/H/2014/home-canning.pdf?la=en
- World Health Organization. Foodborne disease outbreaks: guidelines for investigation and control. Geneva: World Health Organization; 2008. Available from: who.int/foodsafety/publications/foodborne disease/outbreak guidelines.pdf

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.¹

Date of implementation	Description and Change
August 8, 2021	Updated references

Appendix F: Full Revision History

Introduction

Date of Implementation	Description and Change
October 2, 2017	 The title of the Guide has been changed from Public Health Inspector's Guide to the Principles and Practices of Environmental Microbiology to Public Health Inspector's Guide to Environmental Microbiology Laboratory Testing. Removed Contact Directory. Refer clients to the PHO Laboratory Customer Service Centre. Added a note stating PHO Laboratory does not analyze floatation tank/isopod and foot throne samples (no provincial guidelines for microbiological testing or standards available).

Food Sample Analysis

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013:
	 Added links to test directory - test information sheets for sample details (e.g., sample requirements (container, requisition, volume), sample handling and transport, test information, testing frequency, turnaround time and reporting) Updated food bacteriology requisition section to reflect updated requisition (1722- 44 (03/2014)) Updated food-laboratory rejection criteria to reflect acceptance protocol

Environmental Swabs

Date of Implementation	Description and Change
October 2, 2017	 Added links to test directory - test information sheets for sample details (e.g., sample requirements (container, requisition, volume), sample handling and transport, test information, testing frequency, turnaround time and reporting) Updated swab-laboratory rejection criteria to reflect acceptance protocols and added interpretation for swabs.
May 6, 2019	Environmental swabs submitted for R.R.O. 1990, Reg. 562: Food Premises, Section 80 will no longer be accepted for testing, as Section 80 was revoked on July 1, 2018. All applicable areas have been updated to reflect this change.
August 8, 2021	Updated the swab requisition details from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition.

Water Sample Analysis

Date of Implementation	Description and Change
October 2, 2017	 Changes from 2013: Added links to test directory (test information sheets) for sample details (e.g., sample requirements (container, requisition, volume), sample handling and transport, test information, testing frequency, turnaround time and reporting) Removed "Plumbing in the following circumstances" passage from the <i>Drinking Water that fall outside Ontario Drinking Water Regulations</i> chart Updated water-laboratory rejection criteria to reflect acceptance protocols
May 6, 2019	 Changed acceptable level of <i>E. coli</i> for public beaches from 100 CFU/100 mL to 200 CFU/100 mL. Updated the Ministry of the Environment, Conservation and Parks (MECP) name and link. Added table: Water Status Codes and Associated Interpretation Statements.

Date of Implementation	Description and Change
August 8, 2021	 Added the option to submit water samples using <u>Environmental</u> <u>Microbiology Investigation Requisition</u> for etiological agent testing.

Legionella Investigations

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013:
	 Added links to test directory (test information sheets) for sample details (e.g., sample requirements (container, requisition, volume), sample handling and transport, test information, testing frequency, turnaround time and reporting) Added details on environmental sample site recommendations from WHO Added laboratory rejection criteria to reflect acceptance protocols
May 6, 2019	Removed Bacteriological Analysis of Water - Single Sample Requisition for Official Agencies.
August 8, 2021	 Updated the details of the swab requisition from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition.
	Updated the volume requirement for water testing.

Glossary

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013: • Glossary expanded and list of acronyms added.
May 6, 2019	Updated the Ministry of the Environment, Conservation and Parks (MECP) name and link.
August 8, 2021	Added hold time definition.

Acronyms

Date of Implementation	Description and Change
August 8, 2021	Updated MOHLTC to MOH-EH.

Appendix A

Date of Implementation	Description and Change
May 6, 2019	 Removed Bacteriological Analysis of Water - Single Sample Requisition for Official Agencies.

Appendix B

Date of Implementation	Description and Change
August 8, 2021	 Updated the details of the swab requisition from Environmental Bacteriology Swab Test requisition to Environmental Microbiology Investigation Requisition.

Appendix C

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013:
	 Reformatted the chart (there was no change to the content) Added Cyclosporiasis section

Appendix D

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013
	New appendixExpanded the interpretation and added references

Appendix E

Date of implementation	Description and Change
October 2, 2017	Changes from 2013: • Updated references
May 6, 2019	Updated references
August 8, 2021	Updated references

References

Date of Implementation	Description and Change
October 2, 2017	Changes from 2013:
	 Updated references
August 8, 2021	Updated references

References

- Ontario Agency for Health Protection and Promotion (Public Health Ontario). Public health inspector's guide [Internet]. Toronto, ON: Queen's Printer for Ontario; 2019 [cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/public-health-inspectors-guide
- Ontario. Ministry of the Environment, Conservation and Parks. Leading to healthier communities
 and economic prosperity through protecting Ontario's air, land and water [Internet]. Toronto,
 ON: Queen's Printer for Ontario; 2021 [cited 2021 Mar 03]. Available from:
 https://www.ontario.ca/page/ministry-environment-conservation-parks
- Canadian Association for Laboratory Accreditation (CALA). CALA directory of laboratories
 [Internet]. Ottawa, ON: Canadian Association for Laboratory Accreditation; 2021 [cited 2021 Mar
 03]. Available from: http://www.cala.ca/cala_directories.html
- International Organization for Standardization. ISO/IEC 17025:2017: General requirements for the competence of testing and calibration laboratories. Geneva: International Organization for Standardization; 2017.
- United States Food and Drug Administration. Bacteriological analytical manual (BAM) [Internet].
 Silver Spring, MD: Food and Drug Administration; 2021 [cited 2021 Mar 03]. Available from: http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm2006949.htm
- 6. Health Canada. Botulism Reference Service for Canada [Internet]. Ottawa, ON: Government of Canada; 2012 [cited 2021 Mar 03]. Available from: http://www.hc-sc.gc.ca/sr-sr/activ/micro/botulism-eng.php
- Ontario Agency for Health Protection and Promotion (Public Health Ontario). Instructions for
 official agencies submitting water samples to the Public Health Ontario Laboratories [Internet].
 Toronto, ON: Queen's Printer for Ontario; 2012 [cited 2021 Mar 03]. Available from:
 http://www.publichealthontario.ca/en/eRepository/Water_Submission_Instruction_PHOL.pdf
- 8. Ontario Agency for Health Protection and Promotion (Public Health Ontario). Test information index [Internet]. Toronto, ON: Queen's Printer for Ontario; 2021 [cited 2021 Mar 03]. Available from:
 - https://www.publichealthontario.ca/en/ServicesAndTools/LaboratoryServices/Pages/Index.aspx
- 9. Ontario Agency for Health Protection and Promotion (Public Health Ontario). Aerobic colony count (ACC) food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/heterotrophic-plate-count-food

- 10.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Coliform food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/coliform-food
- 11.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Escherichia coli food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/e-coli-food
- 12.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Total gram negative count food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/total-gram-negative-count-food
- 13.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Yeast and mould food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/yeast-mould-food
- 14.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Bacillus cereus food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/bacillus-cereus-food
- 15.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Campylobacter spp. food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/campylobacter-jejuni-food
- 16.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Clostridium perfringens food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/clostridium-perfringens-food
- 17.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Escherichia coli O157:H7 food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2019 Nov 22; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/e-coli-o157-h7-food
- 18.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Listeria monocytogenes food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2019 Nov 22; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/listeria-monocytogenes-food

- 19.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Salmonella spp. food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2019 Nov 22; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/salmonella-food
- 20.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Shigella spp. food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/shigella-food
- 21.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Staphylococcus aureus food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/staph-aureus-food
- 22.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Staphylococcus aureus enterotoxin food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/staph-aureus-toxin-food
- 23.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Vibrio spp. food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/vibrio-food
- 24.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Yersinia food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/yersinia-food
- 25.Ontario Agency for Health Protection and Promotion (Public Health Ontario). pH food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/ph-food
- 26.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Phosphatase food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/phosphatase-food
- 27.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Water activity food [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021

- Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/water-activity-food
- 28.Government of Canada. Canadian Food Inspection Agency [Internet]. Ottawa, ON: Government of Canada; 2021 [modified 2021 Feb 1; cited 2021 Mar 03]. Available from: http://www.inspection.gc.ca/eng/1297964599443/1297965645317
- 29.Ontario. Ministry of Health and Long-Term Care. Ontario's foodborne illness outbreak response protocol (ON-FIORP) 2020: Ontario's guide to a multi-party response to foodborne hazard and illness outbreak investigation [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [cited 2021 Apr 06]. Available from: https://www.health.gov.on.ca/en/pro/programs/publichealth/enviro/docs/fiorp_protocol.pdf
- 30.Health Canada, Health Products and Food Branch (HPFB). Volume 1 Official methods for the microbiological analysis of foods. In: The compendium of analytical methods [Internet]. Evergreen ed. Ottawa, ON: Government of Canada; 1981 [modified 2018 May 09; cited 2021 Mar 03]. Interpretive summary: Health Products and Food Branch (HPFB) standards and guidelines for microbiological safety of food an interpretive summary. Available from: https://www.canada.ca/en/health-canada/services/food-nutrition/research-programs-analytical-methods/analytical-methods/compendium-methods/official-methods-microbiological-analysis-foods-compendium-analytical-methods.html
- 31.Health Canada. Microbial guidelines for ready-to-eat foods a guide for the conveyance industry and Environmental Health Officers (EHO). Ottawa, ON: Her Majesty the Queen in Right of Canada, represented by the Minister of Health; 2013. Available from: publications.gc.ca/collections/collection_2014/sc-hc/H164-167-2013-eng.pdf
- 32.Alberta Health Services. Microbial guidelines for ready-to-eat foods [Internet]. Edmonton, AB: Alberta Health Services; 2011 [modified 2011 May; cited 2017 Jun 06]. Available from: albertahealthservices.ca/assets/wf/eph/wf-eh-microbial-guidelines-for-ready-to-eat-foods.pdf
- 33. Food Standards Australia New Zealand. Compendium of microbiological criteria for food.

 Kingston, Australia: Food Standards Australia New Zealand; 2018. Available from:

 https://www.foodstandards.gov.au/publications/Documents/Compedium%20of%20Microbiological%20Criteria/Compendium_revised-Sep%202018.pdf
- 34.United States Food and Drug Administration. Evaluation and definition of potentially hazardous foods. Silver Spring, MD: United States Food and Drug Administration; 2001. Chapter 3, Factors that influence microbial growth; p. 21-32. Available from: www.fda.gov/downloads/Food/FoodbornelllnessContaminants/UCM545171.pdf
- 35.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Home canning: literature review. Toronto, ON: Queen's Printer for Ontario; 2014. Available from: https://www.publichealthontario.ca/-/media/documents/H/2014/home-canning.pdf?la=en

- 36.Cornell University Department of Food Science. Dairy foods science notes: alkaline phosphatase testing for milk pasteurization [Internet]. Ithaca, NY: Cornell University Department of Food Science; [2007] [cited 2021 Mar 03]. Available from:

 https://foodsafety.foodscience.cornell.edu/files/shared/documents/CU-DFScience-Notes-Milk-Alk-Phosphatase-11-07.pdf
- 37.Manitoba Government. Agriculture: water content and water activity: two factors that affect food safety [Internet]. Winnipeg, MB: Queen's Printer for Manitoba; 2021 [cited 2021 Mar 18]. Available from: https://www.gov.mb.ca/agriculture/food-safety/at-the-food-processor/print,water-content-water-activity.html
- 38.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Drinking water testing official agencies [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2019 Nov 22; cited 2021 Mar 03]. Available from:

 https://www.publichealthontario.ca/en/laboratory-services/test-information-index/drinking-water-testing-official-agencies
- 39.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Drinking water testing private citizen [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from:

 https://www.publichealthontario.ca/en/laboratory-services/test-information-index/drinking-water-testing-private-citizen
- 40.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Ice treated [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/ice-treated
- 41.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Public beach water [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/public-beach-water
- 42.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Recreational water facilities and public pools spas [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from:

 https://www.publichealthontario.ca/en/laboratory-services/test-information-index/recreational-water-facilities-pools-spas
- 43.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Sewage [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Jul 20; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/sewage

44.Ontario. Ministry of the Environment, Laboratory Services Branch. Practices for the collection and handling of drinking water samples [Internet]. Version 2.0. Toronto, ON: Queen's Printer for Ontario; 2009 [cited 2021 Mar 03]. Available from:

https://files.ontario.ca/practices_for_the_collection_and_handling_of_drinking_water_samples.pdf

45. Safe Drinking Water Act, SO 2002, c 32. Available from:

https://www.ontario.ca/laws/statute/02s32

46. Drinking Water Systems, O Reg 170/03. Available from:

https://www.ontario.ca/laws/regulation/030170

47. Small Drinking Water Systems, O Reg 319/08. Available from:

https://www.ontario.ca/laws/regulation/080319

48. Health Protection and Promotion Act, RSO 1990, c H.7. Available from:

https://www.ontario.ca/laws/statute/90h07

49. Camps in Unorganized Territory, RRO 1990, Reg 554. Available from:

https://www.ontario.ca/laws/regulation/900554

50.Ontario. Ministry of the Environment, Laboratory Services Branch. Protocol of accepted drinking water testing methods [Internet]. Version 2.0. Toronto, ON: Queen's Printer for Ontario; 2010 [cited 2021 Mar 03]. Available from:

https://files.ontario.ca/protocol of accepted drinking water testing methods.pdf

51. Drinking Water Testing Services, O Reg 248/03. Available from:

https://www.ontario.ca/laws/regulation/030248

- 52. Ontario Drinking Water Quality Standards. O Reg 169/03. Available from: http://www.e-laws.gov.on.ca/html/regs/english/elaws_regs_030169_e.htm
- 53.Ontario. Ministry of Health and Long-Term Care. Drinking water systems under the *Health Protection and Promotion Act*: directives: get to know the law (II): working together to safeguard our health [Internet]. Toronto, ON: Queen's Printer for Ontario; 2009 [cited 2021 Mar 05]. Available from: https://collections.ola.org/mon/23005/292122.pdf
- 54.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Legionella spp. water [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2020 Dec 09; cited 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/legionella-water
- 55.Ontario Agency for Health Protection and Promotion (Public Health Ontario). Legionella spp. swab [Internet]. Toronto, ON: Queen's Printer for Ontario; 2020 [modified 2019 Nov 22; cited

- 2021 Mar 03]. Available from: https://www.publichealthontario.ca/en/laboratory-services/test-information-index/legionella-swab
- 56. World Health Organization. Legionella and the prevention of legionellosis. Geneva: World Health Organization; 2007. Available from:

 www.who.int/water sanitation health/emerging/legionella.pdf
- 57.Centers for Disease Control and Prevention. *Legionella* environmental assessment form [Internet]. Atlanta, GA: Centers for Disease Control and Prevention; 2019 [cited 2021 Feb 24]. Available from: www.cdc.gov/legionella/downloads/legionella-environmental-assessment.pdf
- 58.Barbaree JM, Gorman GW, Martin WT, Fields BS, Morrill WE. Protocol for sampling environmental sites for legionellae. Appl Environ Microbiol. 1987;53(7):1454-8. Available from: https://www.doi.org/10.1128/AEM.53.7.1454-1458.1987
- 59.Centers for Disease Control and Prevention. Sampling procedure and potential sampling sites: protocol for collecting environmental samples for *Legionella* culture during a cluster or outbreak investigation or when cases of disease may be associated with a facility [Internet]. Atlanta, GA: Centers for Disease Control and Prevention; 2019 [cited 2021 Mar 05]. Available from: https://www.cdc.gov/legionella/downloads/cdc-sampling-procedure.pdf
- 60.Public Health Agency of Canada. Pathogen safety data sheets [Internet]. Ottawa, ON:
 Government of Canada; 2021 [modified 2021 Jan 06; cited 2021 Mar 06]. Available from:
 https://www.canada.ca/en/public-health/services/laboratory-biosafety-biosecurity/pathogen-safety-data-sheets-risk-assessment.html
- 61.Government of Canada. Justice laws website [Internet]. Ottawa, ON: Government of Canada; 2021 [modified 2021 Feb 18; cited 2021 Mar 06]. Available from: http://laws-lois.justice.gc.ca/eng/acts/F-27/
- 62.Ontario. Ministry of Health and Long-Term Care, Population and Public Health Division.

 Operational approaches for recreational water guideline, 2018. Toronto, ON: Queen's Printer for Ontario; 2018. Available from:

 https://www.health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/protocols_guidelines/Operational_Approaches_to_Rec_Water_Guideline_2018_en.pdf
- 63. Canadian Food Inspection Agency. Conducting a hazard analysis [Internet]. Ottawa, ON: Government of Canada; 2019 [modified 2019 Jun 20; cited 2021 Mar 18]. Available from: https://www.inspection.gc.ca/preventive-controls/hazard-analysis/eng/1513283555932/1528205368359
- 64.*Transitional Small Drinking Water Systems,* O Reg 318/08. Available from: http://www.e-laws.gov.on.ca/html/regs/english/elaws_regs_080318_e.htmhttp://aem.asm.org/content/53/7/1454.full.pdf%2Bhtml

- 65. Ashbolt NJ, Grabow WOK, Snozzi M. Indicators of microbial water quality. In: Fewtrell L, Bartram J, editors. Water quality: guidelines, standards and health: assessment of risk and risk management for water-related infectious disease. London: IWA Publishing; 2001. p. 289-316.

 Available from: who.int/water_sanitation_health/dwq/iwachap13.pdf?ua=1
- 66. Municipal Act, 2001, SO 2001, c 25. Available from: https://www.ontario.ca/laws/statute/01m25
- 67. City of Toronto Act, 2006, SO 2006, c 11, Sched. A. Available from: https://www.ontario.ca/laws/statute/06c11
- 68. Provincial Parks and Conservation Reserves Act, 2006, SO 2006, c 12. Available from: ontario.ca/laws/statute/06p12
- 69. Definitions of Words and Expressions Used in the Act, O Reg 171/03. Available from: ontario.ca/laws/regulation/030171
- 70.Ontario. Ministry of Health and Long-Term Care. Recreational water protocol, 2019. Toronto, ON: Queen's Printer for Ontario; 2019. Available from: https://collections.ola.org/mon/2019/05/350606.pdf
- 71. Public Pools, RRO 1990, Reg 565. Available from: ontario.ca/laws/regulation/900565
- 72. Public Spas, O Reg 428/05. Available from: ontario.ca/laws/regulation/050428
- 73.Health Canada, Canadian Food Inspection Agency, Public Health Agency of Canada. Policy on Listeria monocytogenes in ready-to-eat foods [Internet]. Ottawa, ON: Government of Canada; 2011 [cited 2021 Mar 05]. Available from: https://www.canada.ca/content/dam/hc-sc/migration/hc-sc/fn-an/alt_formats/pdf/legislation/pol/policy_listeria_monocytogenes_2011-eng.pdf
- 74. Hong Kong. Centre for Food Safety. Microbiological guidelines for ready-to-eat foods [Internet]. Hong Kong: Centre for Food Safety; 2007 [cited 2021 Mar 05]. Available from: cfs.gov.hk/english/whatsnew/whatsnew act/files/MBGL RTE%20food e.pdf
- 75. World Health Organization. Foodborne disease outbreaks: guidelines for investigation and control. Geneva: World Health Organization; 2008. Available from: who.int/foodsafety/publications/foodborne disease/outbreak guidelines.pdf
- 76.Government of Canada. Cyclosporiasis (Cyclospora) [Internet]. Ottawa, ON: Government of Canada; 2017 [modified 2020 Oct 21; cited 2021 Mar 08]. Available from: https://www.canada.ca/en/public-health/services/diseases/cyclosporiasis-cyclospora.html
- 77. World Health Organization. Guidelines for drinking-water quality. 4th ed. Geneva: World Health Organization; 2017. Chapter 11, Microbial fact sheets; p. 231-306. Available from: apps.who.int/iris/bitstream/10665/254637/1/9789241549950-eng.pdf

- 78.American Public Health Association, American Water Works Association, Water Environment Federation. Standard methods for the examination of water and wastewater. 23rd ed. Washington, DC: American Public Health Association; 2017. Section 9222: Membrane filter technique for members of the coliform group.
- 79.Ontario. Ministry of the Environment, Conservation and Parks. Membrane filtration method using DC agar for the simultaneous detection and enumeration of total coliforms and *Escherichia coli* in drinking water and ground water, Revision 4. Toronto, ON: Queen's Printer for Ontario; 2020.
- 80.Crushell E, Harty S, Sharif F, Bourke B. Enteric campylobacter: purging its secrets? Pediatr Res. 2004;55(1):3-12. Available from: https://doi.org/10.1203/01.pdr.0000099794.06260.71
- 81.Till D, McBride G, Ball A, Taylor K, Pyle E. Large-scale freshwater microbiological study: rationale, results and risks. J Water Health. 2008;6(4):443-60. Available from: https://doi.org/10.2166/wh.2008.071
- 82.Feng P, Weagant SD, Grant MA, Burkhardt W. Enumeration of *Escherichia coli* and the coliform bacteria. In: Tomlinson L (Center for Food Safety and Applied Nutrition, FDA), editor. Bacterial analytical manual. 8th ed. Washington, DC: United States Food and Drug Administration; 2013. Available from: https://www.fda.gov/food/laboratory-methods-food/bam-chapter-4-enumeration-escherichia-coli-and-coliform-bacteria
- 83.Health Canada. Guidelines for Canadian recreational water quality. 3rd edition. Ottawa, ON: Her Majesty the Queen in Right of Canada, represented by the Minister of Health; 2012. Available from: https://healthy-living/guidelines-canadian-recreational-water-quality-third-edition.html?page=9#a411
- 84. Food and Drug Regulations, CRC, c 870. Available from: https://laws.justice.gc.ca/eng/regulations/C.R.C.%2C c. 870/index.html
- 85. Kilian M, Bülow P. Rapid identification of *Enterobacteriaceae*. II. Use of a beta-glucuronidase detecting agar medium (PGUA agar) for the identification of E. coli in primary cultures of urine samples. Acta Pathol Microbiol Scand B. 1979;87(5):271-6.
- 86.Tenaillon O, Skurnik D, Picard B, Denamur E. The population genetics of commensal *Escherichia coli*. Nat Rev Microbiol. 2010;8(3):207-17. Available from: https://doi.org/10.1038/nrmicro2298
- 87.Centre for Food Safety. Microbiological guidelines for food: for ready-to-eat food in general and specific food items [Internet]. Hong Kong: Centre for Food Safety; 2014 [cited 2021 Mar 05]. Available from:
 - cfs.gov.hk/english/food leg/files/food leg Microbiological Guidelines for Food e.pdf
- 88. Dufour AP. Health effects criteria for fresh recreational waters [Internet]. Cincinnati, OH: United States Environmental Protection Agency; 1984 [cited 2020 Mar 05]. Available from:

https://nepis.epa.gov/Exe/ZyNET.exe/300000H7.txt?ZyActionD=ZyDocument&Client=EPA&Index =1981%20Thru%201985&Docs=&Query=&Time=&EndTime=&SearchMethod=1&TocRestrict=n&Toc=&TocEntry=&QField=&QFieldYear=&QFieldMonth=&QFieldDay=&UseQField=&IntQFieldOp=0&ExtQFieldOp=0&XmlQuery=&File=D%3A%5CZYFILES%5CINDEX%20DATA%5C81THRU85%5CTX T%5C00000002%5C300000H7.txt&User=ANONYMOUS&Password=anonymous&SortMethod=h%7C-

&MaximumDocuments=1&FuzzyDegree=0&ImageQuality=r75g8/r75g8/x150y150g16/i425&Display=hpfr&DefSeekPage=x&SearchBack=ZyActionL&Back=ZyActionS&BackDesc=Results%20page&MaximumPages=1&ZyEntry=1

- 89. Wade TJ, Pai N, Eisenberg JN, Colford JM, Jr. Do U.S. Environmental Protection Agency water quality guidelines for recreational waters prevent gastrointestinal illness? A systematic review and meta- analysis. Environ Health Perspect. 2003;111(8):1102-9. Available from: https://doi.org/10.1289/ehp.6241
- 90.Lévesque B, Gauvin D. Microbiological guideline values for recreational bathing in Canada: time for change? Can J Infect Dis Med Microbiol. 2007;18(2):153-7. Available from: https://doi.org/10.1155/2007/180308
- 91.World Health Organization. Guidelines for drinking-water quality. 4th ed. Geneva: World Health Organization; 2017. Chapter 7, Microbial aspects. p. 117-54. Available from: apps.who.int/iris/bitstream/10665/254637/1/9789241549950-eng.pdf
- 92.Neogen Culture Media. Plate Count Agar (Standard Methods) (NCM0010) [Internet]. Lansing, MI: Neogen Corporation; 2019 [cited 2021 Mar 18]. Available from: https://www.neogen.com/globalassets/pim/assets/original/10006/ncm0010_ts_en-us.pdf
- 93.Bartram J, Cotruvo J, Exner M, Fricker C, Glasmacher A, editors. Heterotrophic plate counts and drinking-water safety: the significance of HPCs for water quality and human health. London: IWA Publishing; 2003. Available from: https://www.who.int/water_sanitation_health/dwq/HPCFull.pdf
- 94.United States Environmental Protection Agency. Comprehensive surface water treatment rules quick reference guide: unfiltered systems [Internet]. Washington, DC: United States Environmental Protection Agency; 2010 [cited 2017 Jun 05]. Available from: <a href="https://nepis.epa.gov/Exe/ZyNET.exe/500025PO.TXT?ZyActionD=ZyDocument&Client=EPA&Index=2000+Thru+2005&Docs=&Query=&Time=&EndTime=&SearchMethod=1&TocRestrict=n&Toc=&TocEntry=&QField=&QFieldYear=&QFieldMonth=&QFieldDay=&IntQFieldOp=0&ExtQFieldOp=0&XmlQuery=&File=D%3A%5Czyfiles%5Clndex%20Data%5C00thru05%5CTxt%5C00000019%5C500025PO.txt&User=ANONYMOUS&Password=anonymous&SortMethod=h%7C-&MaximumDocuments=1&FuzzyDegree=0&ImageQuality=r75g8/r75g8/x150y150g16/i425&Display=hpfr&DefSeekPage=x&SearchBack=ZyActionL&Back=ZyActionS&BackDesc=Results%20page&MaximumPages=1&ZyEntry=1&SeekPage=x&ZyPURL

- 95.Food and Drug Regulations, CRC, c 870. s B.12.004. Available from: https://laws.justice.gc.ca/eng/regulations/C.R.C.%2C c. 870/page-64.html#h-572554
- 96. Warburton DW, Dodds KL, Burke R, Johnston MA, Laffey PJ. A review of the microbiological quality of bottled water sold in Canada between 1981 and 1989. Can J Microbiol. 1992;38(1):12-9. Available from: https://doi.org/10.1139/m92-002
- 97.American Public Health Association; American Water Works Association; Water Environment Federation. Standard methods for the examination of water and wastewater. 23rd ed. Washington, DC: American Public Health Association; 2017. Section 9260: Detection of pathogenic bacteria.
- 98. Public Health Agency of Canada. Legionella [Internet]. Ottawa, ON: Government of Canada; 2019 [modified 2019 Nov 26; cited 2021 Mar 05]. Available from: https://www.canada.ca/en/public-health/services/infectious-diseases/legionella.html
- 99. World Health Organization. Legionellosis [Internet]. Geneva: World Health Organization; [cited 2021 Mar 05]. Available from: who.int/mediacentre/factsheets/fs285/en/
- 100. Ontario Agency for Health Protection and Promotion (Public Health Ontario). Legionellosis (legionella, Legionnaires disease) [Internet]. Toronto, ON: Queen's Printer for Ontario; 2019 [cited 2021 Mar 05]. Available from: https://www.publichealthontario.ca/en/diseases-and-conditions/infectious-diseases/respiratory-diseases/legionellosis
- 101. Mena KD, Gerba CP. Risk assessment of *Pseudomonas aeruginosa* in water. Rev Environ Contam Toxicol. 2009;201:71-115. Available from: https://doi.org/10.1007/978-1-4419-0032-6_3
- 102. United States Environmental Protection Agency. Health risks from microbial growth and biofilms in drinking water distribution systems [Internet]. Washington, DC: United States Environmental Protection Agency; 2002 [cited 2021 Feb 24]. Available from: epa.gov/sites/production/files/2015-09/documents/2007_05_18_disinfection_tcr_whitepaper_tcr_biofilms.pdf
- 103. Jacobson JA. Pool-associated *Pseudomonas aeruginosa* dermatitis and other bathing-associated infections. Infect Control. 1985;6(10):398-401. Available from: https://doi.org/10.1017/s0195941700063475
- 104. Botzenhart K, Düring G. Pseudomonas aeruginosa as an opportunistic pathogen. Boston, MA: Springer US; 1993. Chapter 1, Ecology and epidemiology of *Pseudomonas aeruginosa*p. p. 1-18.
- 105. Price D, Ahearn DG. Incidence and persistence of *Pseudomonas aeruginosa* in whirlpools. J Clin Microbiol. 1988;26(9):1650-4. Available from: https://jcm.asm.org/content/26/9/1650.long

- 106. International Organization for Standardization. ISO 16266:2006: Water quality detection and enumeration of *Pseudomonas aeruginosa* method by membrane filtration. Geneva: International Organization for Standardization; 2006.
- 107. Ontario. Ministry of Health and Long-Term Care, Public Health Division, Environmental Health Branch. Response to adverse drinking water quality incidents [Internet]. Toronto, ON: Queen's Printer for Ontario; 2009 [cited 2021 Mar 05]. Available from:

 https://collections.ola.org/mon/23007/293938.pdf

 http://www.boalth.gov.on.ca/on/pro/programs/publishealth/only_standards/docs/guidance/guidanc
 - http://www.health.gov.on.ca/en/pro/programs/publichealth/oph_standards/docs/guidance/guidance_adverse_drinking_water.pdf
- 108. United Kingdom. Health Protection Agency, Centre for Infections. Management of spa pools controlling the risks of infection. London: Health Protection Agency; 2016. Part 2.3.3, Monitoring for infectious agents; p. 56-7. Available from: https://webarchive.nationalarchives.gov.uk/20140714113217/http://www.hpa.org.uk/webc/HPAwebFile/HPAweb C/1200471667554
- 109. Public Health Agency of Canada. Pathogen safety data sheets: infectious substances Salmonella enterica spp. [Internet]. Ottawa, ON: Government of Canada; 2010 [modified 2011 Feb 18; cited 2021 Mar 05]. Available from: phac-aspc.gc.ca/lab-bio/res/psds-ftss/salmonella-ent-eng.php
- 110. International Organization for Standardization. ISO 19250:2010: Water quality detection of *Salmonella* spp. Geneva: International Organization for Standardization; 2010.
- 111. Klapes NA, Vesley D. Rapid assay for in situ identification of coagulase-positive staphylococci recovered by membrane filtration from swimming pool water. Appl Environ Microbiol. 1986;52(3):589-90. Available from: https://doi.org/10.1128/aem.52.3.589-590.1986
- 112. Food Premises, RRO 1990, Reg 562. Available from: https://www.ontario.ca/laws/regulation/900562

Revision History

As an 'evergreen' document the revisions noted below reflect current changes to this Guide. Printed versions of this document should be checked against the <u>Public Health Ontario website</u> to ensure the most recent version is being referenced.

Date of Implementation	Description and Change
August 8, 2021	 Updated references

Public Health Ontario

480 University Avenue, Suite 300 Toronto, Ontario M5G 1V2 647.260.7100 communications@oahpp.ca

publichealthontario.ca